



Overview

Exterior Glove Port Covers may be placed over glove ports with or without gloves attached to facilitate evacuate/fill operations and preserve inert atmosphere. Molded rubber. Includes two clamps. Set of two.

Catalog Number: 5242200

Specifications

- **Weight:** 3.5 lbs
- **Weight metric:** 1.6 kg
- **Applicable Products:** Precise Basic Glove Boxes, Precise Controlled Atmosphere Glove Boxes, Precise HEPA-Filtered Glove Boxes, Protector Combination Glove Boxes, Protector Controlled Atmosphere Glove Boxes, Protector Filtered Glove Boxes, Protector Multi-Hazard Glove Boxes (discontinued), XPert Weigh Boxes

For Product Assistance



[International Customer Service](#)



[Durai Murugan](#)
Regional Sales Manager, India

