



VWR® Roller Bottle

Description

1000mL, Surface Treated, Plug Seal Cap/Vent Cap, Sterile

Purpose

For use in applications where cell attachment is desired or not desired

Materials


Bottle: GPPS (General Polystyrene) Color: Clear
Cap: HDPE (High-density Polyethylene) Color: Green
Membrane: 0.22µm hydrophobic PVDF (only use for Vent Cap)

Features



- Temperature: -20°C ~ +50°C
- Each bottle has silk-printed graduation
- Autoclavability: No
- Working volume: 100-150 mL
- Growth area: 490 cm²
- DNase/RNase free and Non-pyrogenic
- Sterilized by irradiation SAL 10⁻⁶ (ISO11137)
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 clean room environment
- Manufactured under ISO13485 and ISO9001 quality management system

Dimensions (Unit: mm)

Plug seal cap



Vent cap



Printing color: Black



VWR® Roller Bottle

VWR EU Cat. No.	VWR NA Cat. No.	Capacity	Cap Style	Surface Type	Sterile	Packaging Configuration
734-2749	10769-794	1000 mL	Plug seal cap	Surface Treated	Yes	2/Bag, 40/Case
734-2821	75784-530	Cap	Vent cap		Yes	1/Bag, 24/Case



VWR® Roller Bottle

Description

2000mL, Surface Treated, Plug Seal Cap/Vent Cap, Sterile

Purpose

For use in applications where cell attachment is desired or not desired


Materials

Bottle: GPPS (General Polystyrene) Color: Clear
Cap: HDPE (High-density Polyethylene) Color: Green
Membrane: 0.22µm hydrophobic PVDF (only use for Vent Cap)

Features

- Temperature: -20°C ~ +50°C
- Each bottle has silk-printed graduation
- Autoclavability: No
- Working volume: 180-260 mL
- Growth area: 850 cm²
- DNase/RNase free and Non-pyrogenic
- Sterilized by irradiation SAL 10⁻⁶ (ISO11137)
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 clean room environment
- Manufactured under ISO13485 and ISO9001 quality management system

Dimensions (Unit: mm)



VWR® Roller Bottle

VWR EU Cat. No.	VWR NA Cat. No.	Capacity	Cap Style	Surface Type	Sterile	Packaging Configuration
734-2750	10769-796	2000 mL	Plug seal cap	Surface Treated	Yes	20/Bag/Case
	10769-798	2000 mL	Plug seal cap	Surface Treated	Yes	5/Bag, 40/Case
734-2752	10769-800	2000 mL	Vent cap	Surface Treated	Yes	2/Bag, 40/Case
734-3368		2000 mL	Plug seal cap	Surface Treated	Yes	Vertical, 2/Bag, 40/Case
734-3645	76376-716	Cap	Vent cap		Yes	1/Bag, 24/Case



VWR® Roller Bottle

Description

2000mL, Rib Design, Surface Treated, Plug Seal Cap/Vent Cap, Sterile

Purpose

For use in applications where cell attachment is desired or not desired


Materials

Bottle: GPPS (General Polystyrene) Color: Clear
Cap: HDPE (High-density Polyethylene) Color: Green
Membrane: 0.22µm hydrophobic PVDF (only use for Vent Cap)

Features

- Temperature: -20°C ~ +50°C
- Each bottle has silk-printed graduation
- Autoclavability: No
- Working volume: 300-400 mL
- Growth area: 1900 cm²
- DNase/RNase free and Non-pyrogenic
- Sterilized by irradiation SAL 10⁻⁶ (ISO11137)
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 clean room environment
- Manufactured under ISO13485 and ISO9001 quality management system

Dimensions (Unit: mm)



VWR® Roller Bottle

VWR EU Cat. No.	VWR NA Cat. No.	Capacity	Cap Style	Surface Type	Sterile	Packaging Configuration
	10769-802	2000 mL	Plug seal cap	Surface Treated	Yes	5/Bag, 40/Case
734-3369		2000 mL	Plug seal cap	Surface Treated	Yes	Vertical, 2/Bag, 40/Case
734-3645	76376-716	Cap	Vent cap		Yes	1/Bag, 24/Case



VWR® Roller Bottle

Description

5000mL, Surface Treated, Plug Seal Cap/Vent Cap, Sterile

Purpose

For use in applications where cell attachment is desired or not desired


Materials

Bottle: GPPS (General Polystyrene) Color: Clear
Cap: HDPE (High-density Polyethylene) Color: Green
Membrane: 0.22µm hydrophobic PVDF (only use for Vent Cap)


Features

- Temperature: -20°C ~ +50°C
- Each bottle has silk-printed graduation
- Autoclavability: No
- Working volume: 340-510 mL
- Growth area: 1700 cm²
- DNase/RNase free and Non-pyrogenic
- Sterilized by irradiation SAL 10⁻⁶ (ISO11137)
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 clean room environment
- Manufactured under ISO13485 and ISO9001 quality management system


Dimensions (Unit: mm)



Plug seal cap



Vent cap



Printing color: Black

VWR® Roller Bottle

VWR EU Cat. No.	VWR NA Cat. No.	Capacity	Cap Style	Surface Type	Sterile	Packaging Configuration
734-3248	76271-454	5000 mL	Plug seal cap	Surface Treated	Yes	1/Bag, 12/Case
734-3249	76271-456	5000 mL	Vent cap	Surface Treated	Yes	1/Bag, 12/Case
734-2821	75784-530	Cap	Vent cap		Yes	1/Bag, 24/Case



VWR® Roller Bottle

Description

5000mL, Rib Design, Surface Treated, Plug Seal Cap/Vent Cap, Sterile

Purpose

For use in applications where cell attachment is desired or not desired


Materials

Bottle: GPPS (General Polystyrene) Color: Clear
 Cap: HDPE (High-density Polyethylene) Color: Green
 Membrane: 0.22µm hydrophobic PVDF (only use for Vent Cap)

Features

- Temperature: -20°C ~ +50°C
- Each bottle has silk-printed graduation
- Autoclavability: No
- Working volume: 800-1300 mL
- Growth area: 4250 cm²
- DNase/RNase free and Non-pyrogenic
- Sterilized by irradiation SAL 10⁻⁶ (ISO11137)
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 clean room environment
- Manufactured under ISO13485 and ISO9001 quality management system

Dimensions (Unit: mm)



Printing color: Black



VWR® Roller Bottle

VWR EU Cat. No.	VWR NA Cat. No.	Capacity	Cap Style	Surface Type	Sterile	Packaging Configuration
734-3250	76271-458	5000 mL	Plug seal cap	Surface Treated	Yes	1/Bag, 12/Bag/Case
734-3251	76271-460	5000 mL	Vent cap	Surface Treated	Yes	1/Bag, 12/Bag/Case
734-2821	75784-530	Cap	Vent cap		Yes	1/Bag, 24/Case