

Material	20846.326
Material description	Methanol
Grade	GPR RECTAPUR
CAS Number	67-56-1
Molecular formula	H ₃ COH
Molecular mass	32.04

Characteristics	Specifications
Assay	≥ 99.5 %
Appearance	Clear colourless liquid
IR Spectrum	Passes test
Acidity	≤ 0.0002 meq/g
Boiling point	64 - 65 °C
Colouration	≤ 10 APHA
Acetone	≤ 10 ppm
Evaporation residue	≤ 50 ppm
Water	≤ 0.2 %

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Isabelle Habay, Head of Laboratory - Briare
VWR International S.A.S.; Z.I. de Vaugereau; FR-45250
Briare; France

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.