

Product Specification — CO-RE Tip

P/N and Product description

- P/N 235948
- 50µl CO-RE Filter Tip, non-sterile, black, for Microlab STAR Line



Material

- Polypropylen
- Filter tip insert: Polyethylene
- conductive

Packaging & processing

- 96 tips/tray, 5 trays/blister, 5760 tips/case
- Color-code label on tip tray: light orange
- Barcode label contains: unique identification number, Lot No., production date, cavity no.
- Barcode type: Code 128 subset B or C
- Tips usable on CO-RE single channel and 96 Probe Head



Certification / Purity grade

- ISO Class 8 Clean Room production, compliant with ISO 9001, GMP and FDA
- CE IVD conform



- DNase-/RNase-free
- Pyrogen-free (endotoxin-free)
- PCR inhibitor-free

Labeling

Each box of disposable tips provides a label which indicates a barcode, part number (REF) and Lot No.

Measurement

- Tip length: 50.4mm
- Inner diameter tip orifice: 0.4mm

HAMILTON

United Kingdom & Ireland
Tel: +44 (0)121-717-0199

China
Tel: +86-21-6164-6567

France
Tel: +33 (01) 69751616

Italy
Tel: +39-39-689-33-93

Denmark, Norway, Sweden,
Finland
Tel: +46 (0) 8 410 273 73

Germany, Switzerland,
Austria, Benelux
Tel: +49 (089) 552649-0