

Certificate of Conformity

VWR Catalogue Number	See List
Description	Hydrometers, M50
Country of Origin	Manufactured in France
Date of Issue (yyyy-mm-dd)	2022-03-28

270 mm Length, Calibrated at 20°C, 0.001 g/mL Division

Catalogue Number	Range [g/mL]
614-0311	0.800 - 0.850 g/mL
614-0312	0.850 - 0.900 g/mL
614-0313	0.900 - 0.950 g/mL
614-0314	0.950 - 1.000 g/mL
614-0315	1.000 - 1.050 g/mL
614-0316	1.050 - 1.100 g/mL
614-0317	1.100 - 1.150 g/mL
614-0318	1.150 - 1.200 g/mL

Quality System Compliance

Products are manufactured under the **ISO 17025** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Glass
Sterilization	Not Applicable.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	Not Applicable.
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	Not Applicable.
BPA Statement	Not Applicable.
DEHP Statement	Not Applicable.
RoHS	Not Applicable.
REACH Statement	Not Applicable.
Storage Conditions	Ambient temperature.
Shelf Life	No Expiry date.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.