

Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® Spinbar® Stir Bars, Cylindrical
Country of Origin	Manufactured in USA
Date of Issue (yyyy-mm-dd)	2026-03-12

VWR Catalogue Number	Length	Diameter	Coating	Unit
58949-061	12 mm (½")	7.9 mm (5/16 ")	PTFE	Each
				Case of 50
58948-956	19 mm (¾")	7.9 mm (5/16 ")	PTFE	Each
				Case of 50
58949-083	25 mm (1")	7.9 mm (5/16 ")	PTFE	Each
				Case of 50
58948-960	32 mm (1¼")	7.9 mm (5/16 ")	PTFE	Each
				Case of 50
58949-108	38.1 mm (1½")	7.9 mm (5/16 ")	PTFE	Each
				Case of 50
58947-150	44.4 mm (1¾")	7.9 mm (5/16 ")	PTFE	Each
				Case of 50
58948-954	51 mm (2")	7.9 mm (5/16 ")	PTFE	Each
				Case of 50
58949-094	25 mm (1")	9.5 mm (3/8 ")	PTFE	Each
				Case of 50
58948-966	32 mm (1¼")	9.5 mm (3/8 ")	PTFE	Each
				Case of 50
58949-119	38.1 mm (1½")	9.5 mm (3/8 ")	PTFE	Each
				Case of 50
58949-130	51 mm (2")	9.5 mm (3/8 ")	PTFE	Each
				Case of 50
58948-968	54.0 mm (2⅛")	9.5 mm (3/8 ")	PTFE	Each
				Case of 50
58949-152	76 mm (3")	12 mm (1/2")	PTFE	Each
				Case of 50
58948-952	101.6 mm (4")	15.9 mm (5/8 ")	PTFE	Each
				Case of 10
58948-962	127 mm (5")	15.9 mm (5/8 ")	PTFE	Each
				Case of 10

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with [ISO 9001:2015](#) standard.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Molded from Polytetrafluoroethylene (PTFE) over ALNICO 5 Magnet meets US Pharmacopeias Class VI. Compliant with FDA REGULATION 21 CFR 177.1550.
Sterilization	The products can be sterilized in an autoclave at 121°C for a cycle of 20 minutes @15psi.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable
Non-Pyrogenic Statement	Not applicable
Latex Statement	The products do not contain Latex.
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Diethylhexyl phthalate (DEHP) is not included in the manufacture or the formulation of the raw material and are not expected to be present.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr ⁶⁺), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.

REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list is used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.
Storage Conditions	Store at ambient warehouse temperatures, in low humidity conditions, protected from direct sunlight and away from harsh chemicals.
Shelf Life	There is no definite shelf life or expiration for these products. However, products show no degradation when maintained in the recommended storage conditions for up to 5 years.

Disclaimer: VWR'S Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co- manufacturing, or processing activities such as design, packaging, labelling, testing, or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use. This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	March 2026