

Material	20812.320
Material description	Butan-2-ol
Grade	AnalaR NORMAPUR
CAS Number	78-92-2
Molecular formula	$\text{H}_3\text{CCH}_2\text{CH}(\text{OH})\text{CH}_3$
Molecular mass	74.12

Characteristics	Specifications
Assay (on anhydrous substance)	$\geq 99.0 \%$
Acidity	$\leq 0.0005 \text{ meq/g}$
Boiling point	$98.5 - 100.5 \text{ }^\circ\text{C}$
Density (20/4)	$0.805 - 0.809$
Evaporation residue	$\leq 50 \text{ ppm}$
Water	$\leq 0.2 \%$

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Isabelle Habay, Head of Laboratory - Briare
VWR International S.A.S.; Z.I. de Vaugereau; FR-45250
Briare; France

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.