

Certificate of Conformity

VWR Catalogue Number	See List
Description	Filter tips
Country of Origin	Manufactured in USA
Date of Issue (yyyy-mm-dd)	2021-08-18

Bulk, Non-Sterile

Catalogue Number	Tip Style	Volume [μL]	Length [mm]
732-0565	UltraFine, extended	0.1 - 10	38.10
732-0561	UltraFine, graduated	0.1 - 10	31.42
732-0526	Bevel Point	1 - 40	49.66
732-0567	Bevel Point	1 - 50	49.66
732-0521	Bevel Point	1 - 100	49.66
732-0571	Bevel Point	1 - 160	49.66
732-0558	Bevel Point, wide orifice	1 - 200	49.66
732-0539	UltraFine	1 - 200	52.07
732-0535	Bevel Point	1 - 300	52.07
732-0532	UltraFine	100 - 1000	88.90

Racked, Non-Sterile

Catalogue Number	Tip Style	Volume [μL]	Length [mm]
732-0542	UltraFine, extended	0.1 - 10	38.10
732-0515	UltraFine, graduated	0.1 - 10	31.42
732-0529	Bevel Point	1 - 160	49.66
732-0559	Bevel Point, wide orifice	1 - 200	49.66
732-0540	UltraFine	1 - 200	52.07
732-0536	Bevel Point	1 - 300	52.07
732-0533	UltraFine	100 - 1000	88.90

Racked, Sterile

Catalogue Number	Tip Style	Volume [µL]	Length [mm]
732-0543	UltraFine, extended	0.1 - 10	38.10
732-0516	UltraFine, graduated	0.1 - 10	31.42
732-1100	Micro	0.1 - 10	38.10
732-1101	UltraFine, graduated	0.5 - 10	49.66
732-0528	Bevel Point	1 - 40	49.66
732-0525	Bevel Point	1 - 40	49.66
732-0520	Bevel Point	1 - 50	49.66
732-1105	UltraFine, FlexTop, graduated	1 - 100	49.66
732-1103	Standard	1 - 100	49.66
732-1102	Bevel Point, graduated	1 - 100	49.66
732-1107	UltraFine, extended	1 - 100	83.82
732-0523	Bevel Point	1 - 100	49.66
732-0573	Bevel Point	1 - 160	49.66
732-0530	Bevel Point	1 - 160	49.66
732-0541	UltraFine	1 - 200	52.07
732-0610	UltraFine, standard	1 - 200	83.82
732-0544	Bevel Point, wide orifice	1 - 200	49.66
732-0537	Bevel Point	300	52.07
732-0611	UltraFine, extended	1 - 300	71.12
732-1106	UltraFine, FlexTop	100 - 1000	75.95
732-0534	UltraFine	100 - 1000	88.90
613-0752	Wide orifice	100 - 1000	88.90
613-0753	Aerosol barrier	100 - 1200	88.90
732-2894	Sartorius® type	1250	88.90

Individually wrapped, sterile

Catalogue Number	Tip Style	Volume [μL]	Length [mm]
732-0566	UltraFine, extended	0.1 - 10	38.10
732-0524	Bevel Point	1 - 40	49.66
732-0569	Bevel Point	1 - 50	49.66
732-0570	Bevel Point	1 - 100	49.66
732-0572	Bevel Point	1 - 160	49.66
732-0560	Bevel Point, wide orifice	1 - 200	49.66
732-0538	UltraFine	1 - 200	52.07
732-0574	Bevel Point	1 - 300	52.07
732-0531	UltraFine	100 - 1000	88.90

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	100% Virgin Polypropylene. Materials used meet FDA Regulations 21CFR Compliance. Resins are USP Class VI certified.
Sterilization	Product labeled as sterile is gamma or beta irradiated and dose released upon ANSI/AAMI/ISO 11137 recommended practices in effect at the time of validation. Products labeled sterile meets a minimum requirement of 10 ⁻⁶ SAL with a specified dose range of 15.2–32 kGy.
ATP Assay	This product has been tested and is ATP free. The assay detection limit is less than 2x10 ⁻¹² mg/μL.
DNase & RNase Free	This product is free of any detectable DNase/RNase contamination, Human DNA and Protease.

BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Testing is conducted and have been shown to be non-toxic.
Non-Pyrogenic Statement	The acceptance level for product is less than 0.05 EU/mL (LAL Assay).
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenol are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Please see RoHS information below.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2021-07-08 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.
Storage Conditions	Store at room temperature.
Shelf Life	3 years.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.