

Certificate of Conformity

VWR Catalogue Number	See List
Description	Laboratory markers, VWR®
Country of Origin	Manufactured in USA
Date of Issue (yyyy-mm-dd)	2022-03-28

VWR Catalogue Number	Description	Colour
811-0065	Laboratory markers, fine tip	Black
811-0066	Laboratory markers, fine tip	Blue
811-0068	Laboratory markers, fine tip	Green
811-0067	Laboratory markers, fine tip	Red

Quality System Compliance

Products are manufactured under the **ISO 13485:2016** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Cap/Plug:	Celcon, colour concentrate
	Nib:	Acrylic fiber.
	Reservoir:	Bonded acetate fibers.
	Barrel:	Celcon, colour concentrate.
	Ink:	Black ribbon.
Sterilization	Not Applicable.	
ATP Assay	Not Applicable.	
DNase & RNase Free	Not tested.	

BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE)..
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenol are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Product does not contain DEHP.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr ⁶⁺), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2022-01-17 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.
Storage Conditions	General warehouse.
Shelf Life	No shelf life.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.