


Small Alpha® Swab

Alpha Series



Description

The Small Alpha® Swab is an excellent general-purpose swab for cleaning small areas. Its laundered, knitted Alpha® polyester head is extremely clean. A compact handle and rigid internal head core provide support and control. The small tip is ideal for cleaning small surfaces and recesses.

Features

- Cleanroom-laundered, double-layer knitted Alpha® polyester tip over rigid core
- Compact, light green polypropylene handle
- Thermally bonded head
- Soft and nonabrasive
- Excellent chemical resistance

Benefits

- Ultralow nonvolatile residue
- Low particle and fiber generation
- Good sorbency
- No contaminating adhesives
- Broad compatibility with cleaning and sampling solvents

Applications

- Cleaning with solvents such as IPA
- Sampling surfaces for TOC or cleaning validation
- Applying lubricants and other liquids
- Removing excess materials

Products

TX Number	Description	Packaging
TX743B	Small Alpha® Swab	Bag-Within-A-Bag® 500 swabs/bag 5 inner bags of 100 swabs 5 bags/case

*Available in Asia/Pacific



North America
1210 South Park Drive
Kernersville, NC 27284
Tel 336.996.7046
Fax 336.996.2297
www.texwipe.com
info@texwipe.com

Europe/Middle East
Saffierlaan 5
2132 VZ Hoofddorp
The Netherlands
Tel +31 88 1307 400
Fax +31 88 1307 499

Asia/Pacific
50 Tagore Lane
#02-01 Entrepreneur Centre
Singapore 787494
Tel +65 6468 9433
Fax +65 6468 6772

Small Alpha[®] Swab

Physical Characteristics

Head material	Knitted Alpha [®] polyester
Head width	4.0 mm (0.157")
Head thickness	3.6 mm (0.142")
Head length	12.0 mm (0.472")
Handle material	Polypropylene
Handle width	3.0 mm (0.118")
Handle thickness	3.0 mm (0.118")
Handle length	58.3 mm (2.295")
Total swab length	70.3 mm (2.768")
Head bond	Thermal
Handle color	Light green
Design notes	Rigid head core; compact handle

Contamination Characteristics

Property	Measurement	Test Method*
Nonvolatile residue	0.01 mg/swab	TM2: Laboratory Testing For Swabs
IPA extractant		
Ions		
Chloride	0.05 µg/swab	
Sulfate	0.14 µg/swab	
Nitrate	0.08 µg/swab	
Phosphate	0.09 µg/swab	
Fluoride	0.03 µg/swab	
Potassium	0.04 µg/swab	
Calcium	0.01 µg/swab	
Sodium	0.04 µg/swab	
Magnesium	0.04 µg/swab	

*Test procedures available upon request.