

NEW PATENT PENDING TWD SCIENTIFIC DISPOSABLE WEIGHING FUNNELS

Assorted packs and XXL Antistatic
now available



Polypropylene | Antistatic Polypropylene | Borosilicate Glass


avantor™
delivered by **VWR™**

VWR.COM

TWD Scientific, LLC

NEW TWD SCIENTIFIC DISPOSABLE WEIGHING FUNNELS

Ideal for all standard laboratory applications



ASWF-5XXL

(VWR Part Number: **77172-154**)

Antistatic disposable weighing funnel now available in XXL for convenient measurement of larger samples.

ASWF-ASSORT / DPWF-ASSORT / GWF-ASSORT (VWR Part Numbers: **77172-152 / 77172-150 / 77172-148**)

All material types (antistatic polypropylene, standard polypropylene, borosilicate glass) now available to purchase in packs of assorted sizes. Each pack contains 2x extra-small, 4x small, 4x medium, 2x large, and 2x extra-large sized funnels.



VWR Cat. No.	Supplier Cat. No.	Length	Funnel ID	Pouring Neck OD	Max Capacity (mL)	Average Wt. Each (g)	Weight Per Pack (lb)	Fits Volumetric Flask Volume (mL)	#/pk	#pk/cs
14221-812	ASWFXS	38.10mm (1.5 in)	15.24mm (0.6 in)	12.20mm (0.48 in)	1 solid 0.5 liquid	0.65	0.3	96-Well Microtubes	96	12
14221-802	DPWF-PP-XS									
14221-810	ASWF1S	43.18mm (1.7 in)	22.86mm (0.9 in)	7.62mm (0.3 in)	2 solid 2 liquid	1.3	0.15	5, 10, 25, 50	50	22
14221-820	DPWF-PP1-S									
14221-808	ASWF2M	63.50mm (2.5 in)	33.02mm (1.3 in)	10.16mm (0.4 in)	10 solid 7 liquid	3.4	0.4	100, 200, 250	50	8
14221-818	DPWF-PP2-M									
14221-806	ASWF3L	88.90mm (3.5 in)	48.26mm (1.9 in)	17.78mm (0.7 in)	20 solid 10 liquid	7.8	0.35	500, 1000	20	8
14221-816	DPWF-PP3-L									
14221-804	ASWF4XL	114.30mm (4.5 in)	66.04mm (2.6 in)	20.32mm (0.8 in)	60 solid 25 liquid	13.2	0.3	1000, 2000	10	6
14221-814	DPWF-PP4-XL									
77172-152	ASWF-ASSORT	2x extra-small, 4x small, 4x medium, 2x large, 2x extra-large			N/A	N/A	0.5	N/A	14	1
77172-150	DPWF-ASSORT									
77172-154	ASWF5XXL	119.71mm (4.7 in)	XX.XXmm (X.XX in)	28.50mm (1.12 in)	> 270	36.0	1.0	N/A	4	5
14221-800	GWF-BS-XS	42.00mm (1.65 in)	18.00mm (0.7 in)	6.50mm (0.25 in)	1 solid 0.5 liquid	4.3	0.15	96-Well Microtubes	6	14
14221-828	GWF-BS1-S	43.18mm (1.7 in)	22.86mm (0.9 in)	7.62mm (0.3 in)	2 solid 2 liquid	4.6	0.15	5, 10, 25, 50	6	14
14221-826	GWF-BS2-M	63.50mm (2.5 in)	33.02mm (1.3 in)	10.16mm (0.4 in)	10 solid 7 liquid	13.0	0.45	100, 200, 250	6	12
14221-824	GWF-BS3-L	88.90mm (3.5 in)	48.26mm (1.9 in)	17.78mm (0.7 in)	20 solid 10 liquid	27.0	0.75	500, 1000	6	5
14221-822	GWF-BS4-XL	114.30mm (4.5 in)	66.04mm (2.6 in)	20.32mm (0.8 in)	60 solid 25 liquid	49.0	1.3	1000, 2000	6	4
77172-148	GWF-ASSORT	2x extra-small, 4x small, 4x medium, 2x large, 2x extra-large			N/A	N/A	1.5	N/A	14	1

