

Material No.: 0328-07
Revision No.: 0

Product Specification

GMP Manufactured Product

Test	Specification
Appearance (clear, colorless solution)	Passes Test
Normality	0.975 – 1.025
Chloride (Cl)	≤ 5 ppm
ACS – Heavy Metals (as Pb)	≤ 1 ppm
Trace Impurities – Iron (Fe)	≤ 0.5 ppm
Endotoxin Concentration (EU/mL)	≤ 2.0

GMP Manufactured Product

CAUTION: For Manufacturing, processing or repackaging

Made from U.S.P. Purified Water and Sodium Hydroxide pellets, N.F. (Multi-Compendial) which meets BP, Ph.Eur. and JP Chemical Specifications

No Class 1,2,3 or other solvents are used or produced in the manufacturing or purification of the product.

Filtered through a 0.2 micron filter.

Storage Condition: Protect from air to avoid absorption of carbon dioxide.

Packaging Site: Phillipsburg Mfg Ctr & DC



Jamie Ethier
Vice President Global Quality

For questions on this Product Specification please contact Technical Services at 855.282.6867 or +1.610.386.1700

Avantor Performance Materials, LLC

100 Matsonford Rd, Suite 200, Radnor, PA 19087. U.S.A. Phone 610.386.1700

Page 1 of 1