

thermo**scientific**

thermo**scientific**

# Thermo Scientific TSG Series Small Capacity Refrigerators

Innovative refrigeration solutions for  
healthcare and research needs

 **VWR**<sup>™</sup>  
part of avantor

order on  
[VWR.COM](http://VWR.COM)

# Designed with you in mind

The Thermo Scientific TSG Series powered by Phononic<sup>∞</sup>™ solid state refrigeration technology, has been designed to meet the rigorous and demanding needs of healthcare, research and industrial applications where solid performance is required, in support of breastmilk and medication storage, analyzer kits and reagents, and even cell culture media.

The TSG Series brings cutting edge innovation into refrigeration by incorporating solid state technology in place of traditional refrigeration methods—this means that you are not only able to maintain the most demanding storage conditions but you also can enjoy increased storage capacity meaning you can get more done in less space. All while using less energy and putting almost no heat into the work environment.

## Available in two sizes:

- Countertop needs
- Under-counter needs



Under counter refrigerator

Countertop refrigerator



Clear intuitive displays for precise temperature control

## Model specific features:

- Remote alarm contacts on the undercounter refrigerator, supporting requirements of research for monitoring
- Datalogging and download capability on countertop models via USB port to support clinical record management requirements

Warranty for both products is 5 years on the solid state refrigeration system and 2 years parts and labor

## Key features of both models:

- Solid state refrigeration technology
- Whisper quiet operation leading to a more comfortable work environment supporting clinical and healthcare needs
- Increased capacities due to fewer refrigeration components means you can store more in less space
- Energy Star rated performance meaning lower operating costs
- Green refrigerants supporting global initiatives to reduce green house gas emissions
- Full alarm set, including hi/lo temperature, door ajar and power fail to support requirements for medical storage

# Product Specifications

## Thermo Scientific TSG series small capacity undercounter refrigerator

VWR Cat. No.	Settable Range	Capacity cu. ft. (liters)	Electrical	Plug Type*	Doors	Shelves	Defrost	Sound Level	Certifications	Interior Dimensions: (D x W x H)	Exterior Dimensions: (D x W x H)	Shipping Weight lbs. (kg)
76204-742 76204-744	+2°C to +9°C	5.5 (156)	100 - 240v, 50/60hz	NEMA 5-15, NEMA 6-15	1, Solid	3	Auto	43 dBA	UL, CE, REACH, RoHS	19.5 x 20.5 x 23 in. (495x 520x 584 mm)	23.6 x 26 x 32 in. (600 x 660 x 813 mm)	128 (58)
76204-746 76204-748	+2°C to +9°C	5.5 (156)	100 - 240v, 50/60hz	NEMA 5-15, NEMA 6-15	1, Glass	3	Auto	43 dBA	UL, CE, REACH, RoHS	19.5 x 20.5 x 23 in. (495x 520x 584 mm)	23.6 x 26 x 32 in. (600 x 660 x 813 mm)	145 (66)



### Designed for the following industries



Diagnostic Kits & Reagents



Molecular Biology



Vaccines & Pharmaceuticals



Clinical Testing



Industrial testing



Under counter refrigerator



Under counter refrigerator



## Thermo Scientific TSG series small capacity countertop refrigerator

VWR Cat. No.	Settable Range	Capacity cu. ft. (liters)	Electrical	Plug Type*	Doors	Shelves	Defrost	Sound Level	Certifications	Interior Dimensions: (D x W x H)	Exterior Dimensions: (D x W x H)	Shipping Weight lbs. (kg)
76204-740	+2°C to +8°C	1.8 (156)	120v, 60hz	NEMA 5-15	1, Solid	3	Auto	36 dBA	UL, CE,	14.3 x 12.2 x 18.4 in. (363 x 310 x 467 mm)	20.5 x 22.1 x 23 in. (521 x 561 x 584 mm)	83 (38)



### Designed for the following industries



Vaccines & Pharmaceuticals



Clinical Testing



Diagnostic Kits & Reagents

## Customer Installed Accessories

For use with refrigerator size	1.8 cu. ft. (50.9L)	5.5 cu. ft. (156L)
Description	VWR Cat. No.	VWR Cat. No.
Rolling Wire Basket		76045-164
Additional Wire Shelf		76044-688
8.5" Roller Base w/Casters	76266-164	
Pyxis Lock Bracket-Top Mount	76266-166	
Omniceil Flexlock adapter	76266-168	



Smart Tracker wireless monitoring device



Capitol Vial Snappies™ Breast Milk Container System



Thermo Scientific Samco Clicktainer™ Vials and specimen containers