



**EDTA, Disodium Salt,
Dihydrate, Crystal**
BAKER ANALYZED ACS Reagent

Product No. 8993
Specifications current as of:
Jun 8 2011

TEST	SPECIFICATION
Meets A.C.S. Specifications	
Meets Reagent Specifications for testing USP/NF monographs	
Assay (C ₁₀ H ₁₄ N ₂ Na ₂ O ₈ ·2H ₂ O)	99.0 - 101.0 %
Insoluble Matter	max. 0.005 %
pH of 5% Solution at 25°C	4.0 - 6.0
Nitrilotriacetic Acid (by polarography)	max. 0.1 %
Heavy Metals (as Pb)	max. 0.005 %
Iron (Fe)	max. 0.01 %

For Laboratory, Research or Manufacturing Use