

MSDS

HEMATOXYLIN STAIN, HARRIS FORMULA (MODIFIED)

Material Safety Data Sheet

Section 1: Chemical Product and Company Identification

Catalog Number: 3530	
Product Identity: HEMATOXYLIN STAIN, HARRIS FORMULA (MODIFIED)	
Manufacturer's Name: RICCA CHEMICAL COMPANY LLC	Emergency Contact(24 hr) -- CHEMTREC® Domestic: 800-424-9300 International: 703-527-3887
CAGE Code: 0V553	
Address: 448 West Fork Dr Arlington, TX 76012	Telephone Number For Information: 817-461-5601
Date Prepared: 7/29/99	Revision: 2 Last Revised: 09/13/2001 Date Printed: 11/01/2006 8:49:11 am

Section 2. Composition/Information on Ingredients

Component	CAS Registry #	Concentration	ACGIH TLV	OSHA PEL
Hematoxylin	517-28-2	0.4 - 0.6	Not Available Not Available	Not Available Not Available
Ethyl Alcohol (Ethanol)	64-17-5	4 - 5	1000 ppm 1880 mg/m3	1000 ppm 1900 mg/m3
Methanol (Methyl Alcohol)	67-56-1	0.1 - 0.3	200 ppm 262 mg/m3	200 ppm 260 mg/m3
Water, Deionized	7732-18-5	Balance	Not Available Not Available	Not Available Not Available
Potassium Iodate	7758-05-6	< 0.1	Not Available Not Available	Not Available Not Available
Aluminum Ammonium Sulfate Dodecahydrate	7784-26-1	5 - 7	Not Available 2 mg/m3	Not Available Not Available

Section 3: Hazard Identification

Emergency Overview: May cause irritation to the eyes. May be harmful if swallowed. Contact may result in allergic reaction. Wash areas of contact with plenty of water. If ingested, dilute with water, induce vomiting then call a physician.



MSDS

HEMATOXYLIN STAIN, HARRIS FORMULA (MODIFIED)

Target Organs: eyes, skin, respiratory tract, gastrointestinal system, central nervous system, liver, pancreas.

Eye Contact: May cause irritation with burning and stinging with possible damage to the cornea and conjunctiva.

Inhalation: May cause irritation.

Skin Contact: Results in drying and cracking which can lead to secondary infections and dermatitis.

Ingestion: Symptoms can include sleep disorders, hallucinations, distorted perceptions, ataxia, motor function changes, convulsions and tremors, coma, headaches, pulmonary changes, alteration of gastric secretions.

Chronic Effects/Carcinogenicity: None

IARC - No.

NTP - No.

OSHA - No.

Reproductive Information: Reproductive effects cited in 'Registry of Toxic Effects of Chemical Substances' for Methanol (Methyl Alcohol). Reproductive effects cited in 'Registry of Toxic Effects of Chemical Substances' for Ethyl Alcohol (Ethanol).

Teratology (Birth Defect) Information: Mutation data cited in 'Registry of Toxic Effects of Chemical Substances' for Methanol (Methyl Alcohol). Mutation data cited in 'Registry of Toxic Effects of Chemical Substances' for Ethyl Alcohol (Ethanol).

Section 4: First Aid Measures - In all cases, seek qualified evaluation.

Eye Contact: Irrigate immediately with large quantity of water for at least 15 minutes. Call a physician if irritation develops.

Inhalation: Remove to fresh air. Give artificial respiration if necessary. If breathing is difficult, give oxygen.

Skin Contact: Wash areas of contact with soap and water for at least 15 minutes. Call a physician if irritation develops.

Ingestion: Dilute immediately with water or milk. Induce vomiting. Call a physician.

Section 5: Fire Fighting Measures

Flash Point: Not Available.

Method Used: Not Available.

LFL: Not Available.

UFL: Not Available.

Extinguishing Media: Use any means suitable for extinguishing surrounding fire.

Fire & Explosion Hazards: Not considered to be a fire or explosion hazard.

Fire Fighting Instructions: Poisonous gases are produced in fire. Continue to cool containers with water after fire is extinguished. For larger fires, use unmanned hose apparatus, if possible. Consider down wind conditions. Do not release runoff from fire-fighting measures to sewers or waterways.

Fire Fighting Equipment: Wear full protective clothing and NIOSH-approved self-contained breathing apparatus with full facepiece operated in the pressure demand or other positive pressure mode.

Section 6: Accidental Release Measures

Absorb with suitable material and dispose of in accordance with local regulations.

Section 7. Handling and Storage

As with all chemicals, wash hands thoroughly after handling. Avoid contact with eyes and skin. Protect from freezing and physical damage. Store at controlled room temperature, 15 - 30°C.

Safety Storage Code: General

Section 8: Exposure Control/Personal Protection

Engineering Controls: No specific controls are needed. Normal room ventilation is adequate.

Respiratory Protection: Normal room ventilation is adequate.

Skin Protection: Chemical resistant gloves.

Eye Protection: Safety glasses or goggles.



For RICCA, SpectroPure, Red Bird, and Solutions Plus Brands
Emergency Contact (24 hr) -- CHEMTREC®
Domestic: 800-424-9300
International: 703-527-3887

MSDS

HEMATOXYLIN STAIN, HARRIS FORMULA (MODIFIED)

Section 9: Physical and Chemical Properties

Appearance: Dark purple liquid
Odor: very mild alcohol odor
Solubility in Water: Infinite
Specific Gravity: Approximately 1

pH: 2 - 3
Boiling Point(°C): Approximately 98
Melting Point(°C): Approximately -5
Vapor Pressure: Not Applicable.

Section 10: Stability and Reactivity

Chemical Stability: Stable under normal conditions of use and storage.
Incompatibility: Oxidizers, Platinum, Sodium, Potassium Dioxide, Bromine Pentafluoride, acetyl Bromide, Acetyl Chloride.
Hazardous Decomposition Products: Burning may produce oxides of Carbon, Nitrogen and Sulfur.
Hazardous Polymerization: Will not occur.

Section 11. Toxicological Information

LD50, Oral, Rat: 5628 mg/kg (Methanol), details of toxic effects not reported other than lethal dose value; 7060 mg/kg (Ethanol), respiration changes noted.

Section 12. Ecological Information

Ecotoxicological Information: Ethanol has moderate toxicity to aquatic life.
Chemical Fate Information: This material is not expected to significantly bioaccumulate. Ethanol is slightly persistent in water, with a half-life of between 2 to 20 days.

Section 13. Disposal Considerations

Absorb with suitable inert material (vermiculite, dry sand, earth) and place in a chemical waste container for proper disposal in an approved waste disposal facility for incineration in a chemical incinerator equipped with scrubber and afterburner. Always dispose of in accordance with local, state and federal regulations.

Section 14. Transport Information

Part Numbers:

This product is not regulated.

Section 15. Regulatory Information (Not meant to be all inclusive - selected regulation represented)

OSHA Status: These items meet the OSHA Hazard Communication Standard (29 CFR 1910.1200) definition of a hazardous material.

TSCA Status: All components of this solution are listed on the TSCA Inventory or are mixtures (hydrates) of items listed on the TSCA Inventory.

Sara Title III:

Section 302 Extremely Hazardous Substances: Not Applicable.

Section 311/312 Hazardous Categories: Acute, Chronic: Yes Fire, Pressure, Reactivity: No

Section 313 Toxic Chemicals: Not Applicable.

California: None Reported.

Pennsylvania: Methanol (Methyl Alcohol) is listed as an Environmental Hazard on the state's Hazardous Substances List. Ethyl Alcohol (Ethanol) is listed as a Basic Hazard on the state's Hazardous Substances List.

RCRA Status: U154

CERCLA Reportable Quantity: Methanol (Methyl Alcohol) - 5,000 pounds.



For RICCA, SpectroPure, Red Bird, and Solutions Plus Brands
Emergency Contact (24 hr) -- CHEMTREC®
Domestic: 800-424-9300
International: 703-527-3887

MSDS

HEMATOXYLIN STAIN, HARRIS FORMULA (MODIFIED)

WHMIS: Not Applicable.

Section 16. Other Information

NFPA Ratings:

Health: 1 Flammability: 0 Reactivity: 0 Special Notice Key:None

HMIS Ratings:

Health: 1 Flammability: 0 Reactivity: 0 Protective Equipment:B (Protective Eyewear, Gloves)

Rev 1, 9-22-2000: Reformatted from WordPerfect® to Microsoft Word®; (Section 1) Revised emergency telephone number to CHEMTREC® 800-424-9300; (Section 8) revised exposure controls; (Section 15) added Florida and Pennsylvania state references.

Rev 2, 10-09-2001: Reformatted to electronic data format.

When handled properly by qualified personnel, the product described herein does not present a significant health or safety hazard. Alteration of its characteristics by concentration, evaporation, addition of other substances, or other means may present hazards not specifically addressed herein and which must be evaluated by the user. The information furnished herein is believed to be accurate and represents the best data currently available to us. No warranty, expressed or implied, is made and RICCA CHEMICAL COMPANY assumes no legal responsibility or liability whatsoever resulting from its use.

