



Material	84364-300
Material description	Sodium peroxide
Grade	ACS
CAS Number	1313-60-6
Molecular formula	Na <sub>2</sub> O <sub>2</sub>
Molecular mass	77,98

Characteristics	Specifications
Assay	Min. 93,0 %
Heavy metals (as Pb)	Max. 0,002 %
Cl (Chloride)	Max. 0,002 %
PO <sub>4</sub> (Phosphate)	Max. 5 ppm
SO <sub>4</sub> (Sulphate)	Max. 0,001 %
Fe (Iron)	Max. 0,005 %

## Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Michel Boisvert, Head of the Laboratory - Lachine  
VWR Anachemia Canada, 255 rue Norman, Lachine, H8R 1A3,  
QC, Canada

