



Material	74980-760
Material description	di-Potassium hydrogen orthophosphate
Grade	ACS
CAS Number	7758-11-4
Molecular formula	K <sub>2</sub> HPO <sub>4</sub>
Molecular mass	174.18

## Characteristics

## Specifications

Assay	Min. 98.0 %
pH (25°C; 5 %)	8.5 --> 9.6
Heavy metals (ICP-OES)	Max. 5 ppm
Insoluble matter	Max. 0.01 %
Loss on drying (105°C)	Max. 1.0 %
Nitrogen compounds (as N)	Max. 0.001 %
Cl (Chloride)	Max. 0.003 %
SO <sub>4</sub> (Sulphate)	Max. 0.005 %
Fe (Iron)	Max. 0.001 %
Na (Sodium)	Max. 0.05 %

## Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Michel Boisvert, Head of the Laboratory - Lachine  
VWR Anachemia Canada, 255 rue Norman, Lachine, H8R 1A3,  
QC, Canada

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.