



VWR® BASIC PROTECTION INDUSTRIAL BOUFFANT CAPS

As per the manufacturer, the below product meets the following criteria:

North American Catalog No:	76008-524
Manufacture Lot Number:	20180930
Description:	VWR Clip Cap Industrial Red 24in Cs1000
Date of Manufacture:	October 19, 2018
Country of Origin:	China

Quality System Compliance

Products are manufactured under the ISO 9001:2008 & ISO13485:2003 standard. Products are Inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP. Inspection records are reviewed and signed off by qualified personnel for product release.

QC Testing

Inspected according to GB/T 24218.3-2010 method.

Tensile Strength (N/5CM):	GB/T 24218.3-2010 standard	MD = 12.50	CD = 7.20
Elongation Rate (%):	GB/T 24218.3-2010 standard	MD = 100.30	CD = 95.60

Product Specifications

Material:	PP 9g/m ² Red
BSE/TSE Statement:	The above product contains no ingredients of animal origin and no material derived from, or exposed to animals affected by or under quarantine for Transmitting, Animal Spongiform Encephalopathy (TSE)/Bovine Spongiform Encephalopathy (BSE).
Non-Pyrogenic:	N/A
DNase & RNase Free:	N/A
Cytotoxicity:	N/A
Latex-Free Statement:	The above product is latex free.

Signed:



Jamie Ethier
VP Global Quality
VWR, Part of Avantor

Date: January 13, 2020