

## FT902 Immersion Cooler

### with immersion probe for rapid cooling of liquids

JULABO immersion coolers are primarily employed for countercooling open bath circulators to ambient temperature or below as well as for rapidly cooling liquids. JULABO immersion coolers have a small footprint and represent an economic solution preventing tap water usage for cooling. Ease of operation make this model series the ideal instruments for a variety of cooling applications.

### Your advantages

- Rapid cooling of liquids down to low temperatures
- Counter-cooling of liquids in combination with heating circulators
- Environmentally friendly by saving precious tap water
- User-friendly operation and handling
- Compact design, small footprint
- Dry ice substitution
- Integrated temperature control and display
- With external Pt100 sensor (200 x 6 mm dia., stainless steel)



Made in Germany

### Technical Data

Order No.	9650892					
Model series	FT Series					
Category	Immersion Coolers					
Working temperature range (°C)	-90 ... +30					
Temperature stability (°C)	±1					
Temperature Display	LED					
Cooling capacity (Medium Ethanol)	°C	20	10	-20	-40	-80
	kW	0.3	0.27	0.24	0.2	0.07
Refrigerant	R23, R404A					
Ambient temperature	5...35 °C					
Dimensions W x L x H (inch)	15 x 21.7 x 23.6					
Weight (LBS)	110					
Cooling of compressor	Air					
Power requirement V / Hz / A	115/60/7					
Available voltage versions	230 V / 50-60 Hz					
	115 V / 60 Hz					

**JULABO Contact**

JULABO USA Inc.  
884 Marcon Boulevard, Allentown, PA 18109  
Phone: 610-231-0250  
Phone: 800-458-5226 (Toll-Free)  
  
info.us@julabo.com

**JULABO Services**

Product finder, Accessory search, 1 PLUS Warranty,  
Catalog download, Callback service, Operating manuals,  
Safety data sheets, Software and more  
  
[www.julabo.com](http://www.julabo.com)

Technical changes without prior notification. Images may deviate from the original.