


Product Specification

Catalog Number:	39199	
Product Name:	n-Hexane, HPLC Grade, 95+%	
Chemical Abstract No:	110-54-2	
EINECS:	203-777-6	
TSCA:	Yes	
Molecular formula:	H ₃ C (CH ₂) ₄ CH ₃	Structure: 
Formula weight:	86.18	

Typical Data			
Appearance	Clear, colorless liquid	Density	0.659
B.P.	69 °C	Flash point	-23 °C
Viscosity	0.30mPa.s (at 25 oC)	Refractive index	1.3750

Specifications			
Assay(n-Hexane):	> 95.0 %	Water	0.03%
Evaporation residue	5 ppm max		
Filtered to 0.2 micron			

Storage	
Store in tightly sealed containers in a cool, dry location away from direct sunlight.	