

# Transferpette® S

## A perfect fit for every hand!

**BRAND. For lab. For life.®**



Precision workmanship, high grade materials, "Made in Germany" quality – it's what you see on first glance and appreciate every time you use the Transferpette S. The finer points of the design will be evident when the instrument is in your hand. Experience how simple, easy and efficient pipetting can be with the Transferpette S. Single- or multi-channel, there's a Transferpette S pipette for your application.

- **Slim hand grip**—Provides a comfortable and secure fit in any hand
- **One-handed operation**—Perform all functions (set volume, pipette, tip ejection) without changing your grip
- **Volume-change protection**—Simple to operate and clearly shows when the volume setting is locked
- **Slim shaft**—Narrow tubes? With the slim shaft, easily pipette into narrow tubes such as 15 mL conical tubes
- **Easy cleaning**—The shaft tucks inside the handgrip for smooth surfaces, making cleaning easier
- **Low forces**—Short pipetting stroke, easy mounting and tip ejection make pipetting a breeze
- **Universal tip cone** fits tips from most quality manufacturers
- **Completely autoclavable** at 121 °C
- **Color-coded pipette** accents make pipette and tip selection easy
- **Easy mounting and ejection** of tips with O-rings on multi-channel pipettes
- **3-year warranty** for single-channel pipettes

### Transferpette® S pipettes, adjustable volume (fixed volume also available)

Capacity, $\mu$ l (color-coded)	Description	BRANDTECH Cat. No.	VWR Cat. No.
0.1-1	D-1	705868	76383-068
0.1-2.5	D-2.5	705869	76383-086
0.5-10	D-10	705870	76383-084
2-20	D-20	705872	76383-082
5-50	D-50	705873	76383-080
10-100	D-100	705874	76383-078
20-200	D-200	705878	76383-076
100-1000	D-1000	705880	76383-074
500-5000	D-5000	705882	76383-072
1000-10000	D-10000	705884	76383-070

### Transferpette S-8 and -12 multi-channel pipettes

Capacity, $\mu$ l (color-coded)	Description	8-channel		12-channel	
		BRANDTECH Cat. No.	VWR Cat. No.	BRANDTECH Cat. No.	VWR Cat. No.
0.5-10	M8-10	705900	76383-100	705920	76383-110
5-50	M8-50	705906	76383-108	705926	76383-118
10-100	M8-100	705908	76383-106	705928	76383-116
20-200	M8-200	705910	76383-104	705930	76383-114
30-300	M8-300	705912	76383-102	705932	76383-112

### Pipette sets and packages, Transferpette S, adjustable volume, includes tips and shelf/rack (w/sets), benchtop rack (w/packages)

Type	Description	Cat. No.	VWR Cat. No.
Micro	D-2.5, D-10, D-100	705890	76383-088
Mid	D-20, D-200, D-1000	705891	76383-098
Macro	D-1000, D-5000, D-10000	705892	76383-096
Standard	D-10, D-100, D-1000	705893	76383-094
Pipette package 1	D-2.5, D-10, D-20, D-200, D-1000	705896	76383-092
Pipette package 2	D-10, D-100, D-200, D-1000, D-5000	705897	76383-090

