

# Certificate of Quality



<b>Part Number</b>	<b>MFLX40726-70</b>
<b>Lot Number</b>	<b>1374323</b>
<b>Part Description</b>	<b>MF XCNNCTRS PP 1/16 20PK CLNRM</b>
<b>Country of Origin</b>	<b>USA</b>

## Inspection & Testing

- It is the responsibility of the user to conduct tests that are deemed necessary to determine the suitability of this product for any particular use.

## Product Specifications & Compliance

- This material meets the FDA requirements outlined in the Code of Federal Regulations 21 CFR 177.1520(a)(3)(i) and (c)3.1a
- According to our information, all other ingredients used in the formulation meet their respective FDA regulations and 21 CFR 177.1520(b). Specifically, this product meets the FDA requirements for food contact, except for cooking, as listed under conditions of use
- This product has been certified as a USP Class VI Plastic (USP, Biological Reactivity Tests, In Vivo). Physicochemical (USP) and Endotoxin (USP) testing have also been performed.
- Tallow/BSE/TSE We do not use tallow or any other animal derived materials in the manufacture or formulation of this product. This product is Animal Derivative Free.
- This product is compliant to RoHs3 Directive (EU) 2015/863 / (Directive 2011-95/EC RoHs2 amended).
- Shelf life is 5 years from the date of manufacture.
- This product was packaged in an ISO Class 7 cleanroom.
- This product is not approved or intended for, and should not be used for medical, clinical, surgical or other patient-oriented applications.
- It is the responsibility of the user to conduct tests that are deemed necessary to determine the suitability of this product for any particular use.

## Storage

- Storage in original packaging, at ambient conditions and away from direct sunlight, is recommended.

## Certificate Approval

- This certificate was issued and approved by the Avantor Fluid Handling Quality and Regulatory Departments and is valid without signature.