

## RNase Away Certificate of Conformity

Lot Number:23310936

Part Number:7003

Item Description: RNASE AWAY, 1 LITER 5/CS

Release Date: 08/14/2023

Molecular BioProducts Inc., a part of Thermo Fisher Scientific Inc. hereby certifies that this product is effective at removing RNase contamination from common laboratory surfaces. Each lot of RNase AWAY® is tested by an independent laboratory to ensure efficacy and consistency.

### Test Summary

Five pieces of laboratory glassware were intentionally contaminated with RNase. Select pieces were treated with RNase AWAY®, after which the pieces were subjected to a combination of DI water rinses and wiping dry. The pieces were extracted with DEPC water, and the extracts were then added to an RNA standard and incubated. Following incubation, the samples were analyzed by gel electrophoresis and evaluated for degradation of the RNA.

Test results confirm that, when used as directed this lot of RNase AWAY® is effective at removing RNase contamination from laboratory surfaces.

For expiration date, please refer to the label on the product bottle.

Certified by: \_\_\_\_\_



12/13/2023

Ana Ibarria  
Interim Quality Manager

