



Material	BDH1156-4LP
Material description	Reagent Alcohol
Grade	ACS

Characteristics	Specifications
Assay (Ethanol + methanol)	94.0 --> 96.0 %
Assay (2-Propanol)	4.0 --> 6.0 %
Appearance	Clear colourless liquid
Color	Max. 10 APHA
Residue after evaporation	Max. 0.001 %
Water	Max. 0.5 %

## Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Michel Boisvert, Head of the Laboratory - Lachine  
 VWR Anachemia Canada, 255 rue Norman, Lachine, H8R 1A3,  
 QC, Canada

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.