


PDS No. 661960/ 661975	<b>PRODUCT DATA SHEET</b>			Page 1 of 1
Revision 01	Cell Culture Flask, 650 ml, Advanced TC™			 greiner bio-one
	Greiner Item-No. 661960 / 661975			
Valid for Item-No.:	<b>661960</b> (sterile)	<b>661975</b> (sterile)		

1.	Description / Specification	
1.1	Description	Cell Culture Flask, 650 ml, high flask design, canted neck, printed graduation (scale 25-200 ml), writing area, sterile 661960: Standard Cell culture Flasks, Advanced TC™ surface 661975: Filter Cap Cell culture Flasks, Advanced TC™ surface
1.2	Dimensions	Flask: 225.3 x 111.0 x 46.3 mm (length, width, height) 661975: pore size of filter membrane: 0.2 µm
1.3	Volume	Total volume: 650 ml Working volume: 20-85 ml Growth area: 175 cm <sup>2</sup>
1.4	Material / Resin	Flask: PS (Polystyrene), free of heavy metal Cap: HDPE (High Density Polyethylene), free of heavy metal Filter: PET (Polyethylene Terephthalat), PTFE (Polytetrafluorethane), free of heavy metal
1.5	Colour	Flask: clear; print: graduation and writing area white Cap: blue Filter: white
1.6	Sterilization	SAL 10 <sup>-3</sup>
1.7	Quality Control	- <u>Raw Material-Control</u> : physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	- For single use only - Expiry date and Lot-No. printed on bottom of flask

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	+4°C to +37 °C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: <a href="http://www.gbo.com/bioscience">www.gbo.com/bioscience</a> →Products →Literature →Technical Information →Chemical Resistance of Resins
2.6	Shelf life	2 years after month of production
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	4
3.2	Pieces / Box	40
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	Certificate of Quality

4.	Other Information	
		-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision	Date	Date	Date	
-	11 January 2013	14 January 2012	14 January 2012	
Date	Name	Name	Name	
-	S. Kaelberer	Dr. T. Schreiber	A. Schulz	