



a xylem brand

YSI sells this product under agreement with the original manufacturer. YSI assumes no legal liability or responsibility for the accuracy, completeness or usefulness of any information set forth in this MSDS.

SAFETY DATA SHEET ALUMINIUM No. 1 TABLETS

1. Identification

Product identifier

Product name ALUMINIUM No. 1 TABLETS

Product number YAP166, YPM166

Internal identification X-166

Recommended use of the chemical and restrictions on use

Application Testing water

Details of the supplier of the safety data sheet

Supplier PALINTEST USA
1455 JAMIKE AVENUE
ERLANGER
KENTUCKY
41018
USA
859 341 7423
TOLL FREE 800 835 9629
859 341 2106
info@palintestusa.com

DISTRIBUTOR: YSI, INC
ADDRESS 1700/1725 BRANNUM LANE
YELLOW SPRINGS,
OHIO
45387
USA
TELEPHONE (937) 767-7241
E-MAIL MSDSINFO@XYLEMINC.COM

Emergency telephone number

Emergency telephone +44 (0)207 858 1228 (24hr)

National emergency telephone number +44 (0)207 858 1228 (24hr)

Emergency telephone number: Chemtrec® :
US & Canada 800-424-9300
(International) +1 (703) 527-3887

2. Hazard(s) identification

Classification of the substance or mixture

Physical hazards Not Classified

Health hazards Skin Corr. 1B - H314 Eye Dam. 1 - H318 Repr. 1B - H360FD

Environmental hazards Not Classified

Label elements

Pictogram



Signal word Danger

Hazard statements H314 Causes severe skin burns and eye damage.
H360FD May damage fertility. May damage the unborn child.

ALUMINIUM No. 1 TABLETS

Precautionary statements	<p>P201 Obtain special instructions before use.</p> <p>P202 Do not handle until all safety precautions have been read and understood.</p> <p>P260 Do not breathe vapor/spray.</p> <p>P264 Wash contaminated skin thoroughly after handling.</p> <p>P280 Wear protective gloves/protective clothing/eye protection/face protection.</p> <p>P301+P330+P331 If swallowed: Rinse mouth. Do NOT induce vomiting.</p> <p>P303+P361+P353 If on skin (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.</p> <p>P304+P340 If inhaled: Remove person to fresh air and keep comfortable for breathing.</p> <p>P305+P351+P338 If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.</p> <p>P308+P313 If exposed or concerned: Get medical advice/attention.</p> <p>P310 Immediately call a poison center/doctor.</p> <p>P321 Specific treatment (see medical advice on this label).</p> <p>P363 Wash contaminated clothing before reuse.</p> <p>P405 Store locked up.</p> <p>P501 Dispose of contents/container in accordance with national regulations.</p>
---------------------------------	---

Contains POTASSIUM HYDROGENSULPHATE, BORIC ACID

Other hazards

This substance is not classified as PBT or vPvB according to current EU criteria.

3. Composition/information on ingredients

Mixtures

ADIPIC ACID	10-30%
CAS number: 124-04-9	REACH registration number: 01-2119457561-38-XXXX

Classification

Eye Irrit. 2 - H319

BORIC ACID	10-30%
CAS number: 10043-35-3	REACH registration number: 01-2119486683-25-XXXX

Classification

Repr. 1B - H360FD

POTASSIUM HYDROGENSULPHATE	10-30%
CAS number: 7646-93-7	REACH registration number: N/A

Classification

Skin Corr. 1B - H314
 Eye Dam. 1 - H318
 STOT SE 3 - H335

The Full Text for all Hazard Statements are Displayed in Section 16.

Composition comments The Data Shown is in accordance with the latest EC Directives.

4. First-aid measures

ALUMINIUM No. 1 TABLETS

Description of first aid measures

Inhalation	Unlikely route of exposure as the product does not contain volatile substances.
Ingestion	Rinse mouth thoroughly with water. Get medical attention immediately.
Skin Contact	Wash skin thoroughly with soap and water. Get medical attention promptly if symptoms occur after washing.
Eye contact	Rinse immediately with plenty of water. Remove any contact lenses and open eyelids wide apart. Continue to rinse for at least 15 minutes. Get medical attention if irritation persists after washing. Show this Safety Data Sheet to the medical personnel.

Most important symptoms and effects, both acute and delayed

General information	The severity of the symptoms described will vary dependent on the concentration and the length of exposure.
Inhalation	This is unlikely to occur but symptoms similar to those of ingestion may develop.
Ingestion	May cause chemical burns in mouth and throat. May cause stomach pain or vomiting.
Skin contact	May cause serious chemical burns to the skin.
Eye contact	May cause blurred vision and serious eye damage.

Indication of immediate medical attention and special treatment needed

Notes for the doctor	Treat symptomatically.
-----------------------------	------------------------

5. Fire-fighting measures

Extinguishing media

Suitable extinguishing media	Use fire-extinguishing media suitable for the surrounding fire.
-------------------------------------	---

Special hazards arising from the substance or mixture

Hazardous combustion products	Heating may generate the following products: Toxic and corrosive gases or vapors.
--------------------------------------	---

Advice for firefighters

Protective actions during firefighting	No specific firefighting precautions known.
---	---

6. Accidental release measures

Personal precautions, protective equipment and emergency procedures

Personal precautions	Wear protective clothing as described in Section 8 of this safety data sheet.
-----------------------------	---

Environmental precautions

Environmental precautions	Not considered to be a significant hazard due to the small quantities used. The product components are not classified as environmentally hazardous. However, large or frequent spills may have hazardous effects on the environment.
----------------------------------	--

Methods and material for containment and cleaning up

Methods for cleaning up	Absorb in vermiculite, dry sand or earth and place into containers. Wear suitable protective equipment, including gloves, goggles/face shield, respirator, boots, clothing or apron, as appropriate. Collect and place in suitable waste disposal containers and seal securely. For waste disposal, see Section 13.
--------------------------------	---

Reference to other sections	Wear protective clothing as described in Section 8 of this safety data sheet.
------------------------------------	---

ALUMINIUM No. 1 TABLETS

7. Handling and storage

Precautions for safe handling

Usage precautions Good personal hygiene procedures should be implemented.

Conditions for safe storage, including any incompatibilities

Storage precautions Keep separate from food, feedstuffs, fertilizers and other sensitive material. Store in tightly-closed, original container in a dry and cool place.

Specific end uses(s)

Specific end use(s) The identified uses for this product are detailed in Section 1.2.

8. Exposure Controls/personal protection

Control parameters

Occupational exposure limits

ADIPIC ACID

Long-term exposure limit (8-hour TWA): ACGIH 5 mg/m³

BORIC ACID

Long-term exposure limit (8-hour TWA): ACGIH 2 mg/m³

Short-term exposure limit (15-minute): ACGIH 6 mg/m³

A4

ACGIH = American Conference of Governmental Industrial Hygienists.

A4 = Not Classifiable as a Human Carcinogen.

Ingredient comments No exposure limits known for ingredient(s).

ADIPIC ACID (CAS: 124-04-9)

Ingredient comments No exposure limits known for ingredient(s).

Exposure controls

Protective equipment



Eye/face protection

Eyewear complying with an approved standard should be worn if a risk assessment indicates eye contact is possible. The following protection should be worn: Chemical splash goggles.

Hand protection

It is recommended that chemical-resistant, impervious gloves are worn. It is recommended that gloves are made of the following material: Nitrile rubber.

Hygiene measures

No specific hygiene procedures recommended but good personal hygiene practices should always be observed when working with chemical products. Wash promptly if skin becomes contaminated.

Respiratory protection

No specific recommendations. Respiratory protection may be required if excessive airborne contamination occurs.

9. Physical and Chemical Properties

Information on basic physical and chemical properties

Appearance Solid

ALUMINIUM No. 1 TABLETS

Color	White.
Odor	Odorless.
Other information	No data available.

10. Stability and reactivity

Reactivity	There are no known reactivity hazards associated with this product.
Stability	Stable under the prescribed storage conditions.
Possibility of hazardous reactions	No data available.
Conditions to avoid	There are no known conditions that are likely to result in a hazardous situation.
Materials to avoid	No specific material or group of materials is likely to react with the product to produce a hazardous situation.
Hazardous decomposition products	There are no anticipated hazardous decomposition products associated with this material.

11. Toxicological information

Information on toxicological effects

Toxicological effects	This product is toxic.
------------------------------	------------------------

12. Ecological Information

Ecotoxicity	There are no data on the ecotoxicity of this product. The product components are not classified as environmentally hazardous. However, large or frequent spills may have hazardous effects on the environment.
--------------------	--

Toxicity

Toxicity	No data available.
-----------------	--------------------

Persistence and degradability

Persistence and degradability	There are no data on the degradability of this product.
--------------------------------------	---

Bioaccumulative potential

Bio-Accumulative Potential	No data available on bioaccumulation.
-----------------------------------	---------------------------------------

Mobility in soil

Mobility	The product is soluble in water.
-----------------	----------------------------------

Results of PBT and vPvB assessment

Results of PBT and vPvB assessment	This substance is not classified as PBT or vPvB according to current EU criteria.
---	---

Other adverse effects

Other adverse effects	None known.
------------------------------	-------------

13. Disposal considerations

Waste treatment methods

ALUMINIUM No. 1 TABLETS

Disposal methods Dispose of waste to licensed waste disposal site in accordance with the requirements of the local Waste Disposal Authority.

14. Transport information

Road transport notes The product is not classified as packaged as per ADR 2013 section 3.5.1.4

Sea transport notes The product is not classified as packaged as per IMDG 2012 section 3.5.1.4

Air transport notes The product is not subject to the requirements of IATA as packaged based upon "De Minimis Quantities" as per IATA Dangerous Goods Regulations (54th edition) subsection 2.6.10.1

UN Number

UN No. (DOT) 2509

UN No. (IMDG) 2509

UN No. (ICAO) 2509

UN proper shipping name

Proper shipping name (DOT) POTASSIUM HYDROGEN SULPHATE

Proper shipping name (IMDG) POTASSIUM HYDROGEN SULPHATE

Proper shipping name (ICAO) POTASSIUM HYDROGEN SULPHATE

Transport hazard class(es)

IMDG Class 8

IMDG subsidiary risk

ICAO class/division 8

ICAO subsidiary risk

Transport labels



Packing group

DOT pack group II

IMDG packing group II

ICAO packing group II

Environmental hazards

Environmentally Hazardous Substance

No.

Special precautions for user

EmS F-A, S-B

Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code Not applicable.

15. Regulatory information

ALUMINIUM No. 1 TABLETS

US Federal Regulations

CERCLA/Superfund, Hazardous Substances/Reportable Quantities (EPA)

None of the ingredients are listed.

SARA 313 Emission Reporting

None of the ingredients are listed.

CAA Accidental Release Prevention

None of the ingredients are listed.

US State Regulations

California Proposition 65 Carcinogens and Reproductive Toxins

None of the ingredients are listed.

Massachusetts "Right To Know" List

None of the ingredients are listed.

Rhode Island "Right To Know" List

None of the ingredients are listed.

Minnesota "Right To Know" List

None of the ingredients are listed.

New Jersey "Right To Know" List

None of the ingredients are listed.

Pennsylvania "Right To Know" List

None of the ingredients are listed.

Inventories

Canada - DSL/NDSL

All the ingredients are listed or exempt.

US - TSCA

All the ingredients are listed or exempt.

US - TSCA 12(b) Export Notification

None of the ingredients are listed.

16. Other information

Issued by	L. Morgan
Revision date	5/21/2015
Revision	10
Supersedes date	25/07/2014
SDS status	Approved.
Hazard statements in full	H314 Causes severe skin burns and eye damage. H318 Causes serious eye damage. H319 Causes serious eye irritation. H335 May cause respiratory irritation. H360FD May damage fertility. May damage the unborn child.

ALUMINIUM No. 1 TABLETS

This information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process. Such information is, to the best of the company's knowledge and belief, accurate and reliable as of the date indicated. However, no warranty, guarantee or representation is made to its accuracy, reliability or completeness. It is the user's responsibility to satisfy himself as to the suitability of such information for his own particular use.



a xylem brand

YSI sells this product under agreement with the original manufacturer. YSI assumes no legal liability or responsibility for the accuracy, completeness or usefulness of any information set forth in this MSDS.

SAFETY DATA SHEET ALUMINIUM No. 2 TABLETS

1. Identification

Product identifier

Product name ALUMINIUM No. 2 TABLETS

Product number YAP166, YPM166

Internal identification X-167

Recommended use of the chemical and restrictions on use

Application Testing water

Details of the supplier of the safety data sheet

Supplier PALINTEST USA
1455 JAMIKE AVENUE
ERLANGER
KENTUCKY
41018
USA
859 341 7423
TOLL FREE 800 835 9629
859 341 2106
info@palintestusa.com

DISTRIBUTOR: YSI, INC
ADDRESS 1700/1725 BRANNUM LANE
YELLOW SPRINGS,
OHIO
45387
USA
TELEPHONE (937) 767-7241
E-MAIL MSDSINFO@XYLEMINC.COM

Emergency telephone number

Emergency telephone +44 (0)207 858 1228 (24hr)

National emergency telephone number +44 (0)207 858 1228 (24hr)

Emergency telephone number: Chemtrec® :
US & Canada) 800-424-9300
(International) +1 (703) 527-3887

2. Hazard(s) identification

Classification of the substance or mixture

Physical hazards Not Classified

Health hazards Skin Sens. 1 - H317

Environmental hazards Not Classified

Label elements

Pictogram



Signal word Warning

Hazard statements H317 May cause an allergic skin reaction.

ALUMINIUM No. 2 TABLETS

Precautionary statements	<p>P261 Avoid breathing vapor/spray.</p> <p>P272 Contaminated work clothing must not be allowed out of the workplace.</p> <p>P280 Wear protective gloves/protective clothing/eye protection/face protection.</p> <p>P302+P352 If on skin: Wash with plenty of water.</p> <p>P321 Specific treatment (see medical advice on this label).</p> <p>P333+P313 If skin irritation or rash occurs: Get medical advice/attention.</p> <p>P362+P364 Take off contaminated clothing and wash it before reuse.</p> <p>P501 Dispose of contents/container in accordance with national regulations.</p>
---------------------------------	--

Contains HEXAMINE

Other hazards

This substance is not classified as PBT or vPvB according to current EU criteria.

3. Composition/information on ingredients

Mixtures

SODIUM HYDROGEN CARBONATE	60-80%
CAS number: 144-55-8	REACH registration number: 01-2119457606-32-XXXX
Classification Not Classified	
HEXAMINE	30-60%
CAS number: 100-97-0	REACH registration number: 01-2119474895-20-XXXX
Classification Flam. Sol. 2 - H228 Skin Sens. 1 - H317	
POTASSIUM CHLORIDE	5-10%
CAS number: 7447-40-7	REACH registration number: 01-2119539416-36-XXXX
Classification Not Classified	
ERIOCHROME CYANINE	<0.5%
CAS number: 3564-18-9	
Classification Not Classified	

The Full Text for all Hazard Statements are Displayed in Section 16.

Composition comments The Data Shown is in accordance with the latest EC Directives.

4. First-aid measures

Description of first aid measures

General information Treat symptomatically.

ALUMINIUM No. 2 TABLETS

Inhalation	Remove affected person from source of contamination. Move affected person to fresh air at once. Get medical attention if any discomfort continues.
Ingestion	Rinse mouth thoroughly with water. Get medical attention if any discomfort continues.
Skin Contact	Remove affected person from source of contamination. Remove contaminated clothing. Wash skin thoroughly with soap and water. Get medical attention if irritation persists after washing.
Eye contact	Remove any contact lenses and open eyelids wide apart. Continue to rinse for at least 15 minutes. Get medical attention promptly if symptoms occur after washing.

Most important symptoms and effects, both acute and delayed

General information	The severity of the symptoms described will vary dependent on the concentration and the length of exposure.
Inhalation	Irritation of nose, throat and airway.
Ingestion	May cause discomfort if swallowed.
Skin contact	Skin irritation.
Eye contact	May cause severe eye irritation.

Indication of immediate medical attention and special treatment needed

Notes for the doctor	Treat symptomatically.
-----------------------------	------------------------

5. Fire-fighting measures

Extinguishing media

Suitable extinguishing media Use fire-extinguishing media suitable for the surrounding fire.

Special hazards arising from the substance or mixture

Hazardous combustion products No known hazardous decomposition products.

Advice for firefighters

Protective actions during firefighting No specific firefighting precautions known.

6. Accidental release measures

Personal precautions, protective equipment and emergency procedures

Personal precautions Wear protective clothing as described in Section 8 of this safety data sheet.

Environmental precautions

Environmental precautions Not considered to be a significant hazard due to the small quantities used. However, large or frequent spills may have hazardous effects on the environment.

Methods and material for containment and cleaning up

Methods for cleaning up Avoid generation and spreading of dust. Collect and place in suitable waste disposal containers and seal securely. Label the containers containing waste and contaminated materials and remove from the area as soon as possible. Flush contaminated area with plenty of water. Wear suitable protective equipment, including gloves, goggles/face shield, respirator, boots, clothing or apron, as appropriate.

Reference to other sections Wear protective clothing as described in Section 8 of this safety data sheet.

7. Handling and storage

ALUMINIUM No. 2 TABLETS

Precautions for safe handling

Usage precautions Good personal hygiene procedures should be implemented. Avoid inhalation of dust and contact with skin and eyes.

Conditions for safe storage, including any incompatibilities

Storage precautions Store in tightly-closed, original container in a dry, cool and well-ventilated place. Keep separate from food, feedstuffs, fertilizers and other sensitive material.

Specific end uses(s)

Specific end use(s) The identified uses for this product are detailed in Section 1.2.

8. Exposure Controls/personal protection

Ingredient comments WEL = Workplace Exposure Limits No exposure limits known for ingredient(s).

Exposure controls

Protective equipment



Appropriate engineering controls Provide adequate ventilation. Avoid inhalation of dust. Observe any occupational exposure limits for the product or ingredients.

Eye/face protection Eyewear complying with an approved standard should be worn if a risk assessment indicates eye contact is possible. The following protection should be worn: Dust-resistant, chemical splash goggles.

Hand protection Wear protective gloves. It is recommended that gloves are made of the following material: Nitrile rubber.

Other skin and body protection Wear appropriate clothing to prevent any possibility of skin contact.

Hygiene measures No specific hygiene procedures recommended but good personal hygiene practices should always be observed when working with chemical products.

Respiratory protection If ventilation is inadequate, suitable respiratory protection must be worn.

9. Physical and Chemical Properties

Information on basic physical and chemical properties

Appearance Solid

Color Beige. Brownish.

Odor Odorless.

Other information No data available.

10. Stability and reactivity

Reactivity There are no known reactivity hazards associated with this product.

Stability Stable under the prescribed storage conditions.

Possibility of hazardous reactions No data available.

ALUMINIUM No. 2 TABLETS

Conditions to avoid	There are no known conditions that are likely to result in a hazardous situation.
Materials to avoid	No specific material or group of materials is likely to react with the product to produce a hazardous situation.
Hazardous decomposition products	No known hazardous decomposition products.

11. Toxicological information

Information on toxicological effects

Toxicological effects No information available.

12. Ecological Information

Ecotoxicity There are no data on the ecotoxicity of this product. The product components are not classified as environmentally hazardous. However, large or frequent spills may have hazardous effects on the environment.

Toxicity

Toxicity No data available.

Persistence and degradability

Persistence and degradability No data available.

Bioaccumulative potential

Bio-Accumulative Potential No data available on bioaccumulation.

Mobility in soil

Mobility The product is soluble in water.

Results of PBT and vPvB assessment

Results of PBT and vPvB assessment This substance is not classified as PBT or vPvB according to current EU criteria.

Other adverse effects

Other adverse effects None known.

13. Disposal considerations

Waste treatment methods

Disposal methods Dispose of waste to licensed waste disposal site in accordance with the requirements of the local Waste Disposal Authority.

14. Transport information

General The product is not covered by international regulations on the transport of dangerous goods (IMDG, IATA, DoT).

UN Number

Not applicable.

UN proper shipping name

Not applicable.

Transport hazard class(es)

ALUMINIUM No. 2 TABLETS

No transport warning sign required.

Packing group

Not applicable.

Environmental hazards

Environmentally Hazardous Substance

No.

Special precautions for user

Not applicable.

Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code Not applicable.

15. Regulatory information

US Federal Regulations

CERCLA/Superfund, Hazardous Substances/Reportable Quantities (EPA)

None of the ingredients are listed.

SARA 313 Emission Reporting

None of the ingredients are listed.

CAA Accidental Release Prevention

None of the ingredients are listed.

US State Regulations

California Proposition 65 Carcinogens and Reproductive Toxins

None of the ingredients are listed.

Massachusetts "Right To Know" List

None of the ingredients are listed.

Rhode Island "Right To Know" List

None of the ingredients are listed.

Minnesota "Right To Know" List

None of the ingredients are listed.

New Jersey "Right To Know" List

None of the ingredients are listed.

Pennsylvania "Right To Know" List

None of the ingredients are listed.

Inventories

Canada - DSL/NDSL

All the ingredients are listed or exempt.

US - TSCA

All the ingredients are listed or exempt.

US - TSCA 12(b) Export Notification

ALUMINIUM No. 2 TABLETS

None of the ingredients are listed.

16. Other information

Issued by	L. Morgan
Revision date	5/21/2015
Revision	5
Supersedes date	25/07/2014
SDS status	Approved.
Hazard statements in full	H317 May cause an allergic skin reaction.

This information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process. Such information is, to the best of the company's knowledge and belief, accurate and reliable as of the date indicated. However, no warranty, guarantee or representation is made to its accuracy, reliability or completeness. It is the user's responsibility to satisfy himself as to the suitability of such information for his own particular use.