



VWR ANALYTICAL

Product Specification

Catalog Number: BDH1158-4LP

Catalog Description: Alcohol 95%

Grade:

Date Current Specification Set: March 12 2012

All batches of this product manufactured after this date are tested to this specification. A new certificate will be published on the web in the unlikely event of specification changes.

Characteristics	Specifications
Appearance:	Clear, colorless liquid
Water:	4.0 – 6.0%
Residue:	Less than 0.001% (0.0010g/122mL)
GC Purity:	
Ethanol	85.0-87.0%
Methanol	4.5-5.0%
Isopropyl Alcohol	4.5-5.5%
Color (APHA):	Less than 10
Specific gravity @ 20°C:	0.806-0.813
Acidity:	Less than 0.0025%

Signed on behalf of VWR:

Neal Fox
Quality Control Supervisor