

2-Propanol  
CMOS  
(iso-propyl alcohol)



Material No.: 9059-09

Revision No.: 0

## Product Specification

Test	Specification
Assay (CH <sub>3</sub> CHOHCH <sub>3</sub> )	≥ 99.5 %
Color (APHA)	≤ 10
Residue after Evaporation	≤ 4 ppm
Solubility in H <sub>2</sub> O	Passes Test
Water (H <sub>2</sub> O)(by Karl Fischer titrn)	≤ 0.05 %
Acidity (μeq/g)	≤ 0.2
Alkalinity (μeq/g)	≤ 0.1
Chloride (Cl)	≤ 0.2 ppm
Phosphate (PO <sub>4</sub> )	≤ 0.3 ppm
Trace Impurities - Aluminum (Al)	≤ 50.0 ppb
Arsenic and Antimony (as As)	≤ 10.0 ppb
Trace Impurities - Barium (Ba)	≤ 20.0 ppb
Trace Impurities - Beryllium (Be)	≤ 100.0 ppb
Trace Impurities - Bismuth (Bi)	≤ 100.0 ppb
Trace Impurities - Boron (B)	≤ 10.0 ppb
Trace Impurities - Cadmium (Cd)	≤ 20.0 ppb
Trace Impurities - Calcium (Ca)	≤ 50.0 ppb
Trace Impurities - Chromium (Cr)	≤ 20.0 ppb
Trace Impurities - Cobalt (Co)	≤ 20 ppb
Trace Impurities - Copper (Cu)	≤ 10.0 ppb
Trace Impurities - Gallium (Ga)	≤ 50 ppb
Trace Impurities - Germanium (Ge)	≤ 50.0 ppb
Trace Impurities - Gold (Au)	≤ 20 ppb
Heavy Metals (as Pb)	≤ 200.0 ppb
Trace Impurities - Iron (Fe)	≤ 50.0 ppb
Trace Impurities - Lead (Pb)	≤ 20.0 ppb
Trace Impurities - Lithium (Li)	≤ 50.0 ppb
Trace Impurities - Magnesium (Mg)	≤ 20 ppb

>>> Continued on page 2 >>>

For questions on this Product Specification please contact Technical Services at 855.282.6867 or +1.610.386.1700

Avantor Performance Materials, LLC

100 Mansford Rd, Suite 200, Radnor, PA 19087. U.S.A. Phone 610.386.1700

2-Propanol  
CMOS  
(iso-propyl alcohol)



Material No.: 9059-09

Test	Specification
Trace Impurities - Manganese (Mn)	≤ 15.0 ppb
Trace Impurities - Molybdenum (Mo)	≤ 100.0 ppb
Trace Impurities - Nickel (Ni)	≤ 10.0 ppb
Trace Impurities - Niobium (Nb)	≤ 100.0 ppb
Trace Impurities - Potassium (K)	≤ 100 ppb
Trace Impurities - Silicon (Si)	≤ 50 ppb
Trace Impurities - Silver (Ag)	≤ 20.0 ppb
Trace Impurities - Sodium (Na)	≤ 100.0 ppb
Trace Impurities - Strontium (Sr)	≤ 20.0 ppb
Trace Impurities - Tantalum (Ta)	≤ 100.0 ppb
Trace Impurities - Thallium (Tl)	≤ 10.0 ppb
Trace Impurities - Tin (Sn)	≤ 100 ppb
Trace Impurities - Titanium (Ti)	≤ 20 ppb
Trace Impurities - Vanadium (V)	≤ 100.0 ppb
Trace Impurities - Zinc (Zn)	≤ 50 ppb
Trace Impurities - Zirconium (Zr)	≤ 100.0 ppb
Particle Count at point of fill - 0.5 µm and greater (Rion KS42AF)	≤ 150 par/ml
Particle Count at point of fill - 1.0 µm and greater (Rion KS42AF)	≤ 25 par/ml

>>> Continued on page 3 >>>

For questions on this Product Specification please contact Technical Services at 855.282.6867 or +1.610.386.1700  
Avantor Performance Materials, LLC  
100 Mansford Rd, Suite 200, Radnor, PA 19087. U.S.A. Phone 610.386.1700

2-Propanol  
CMOS  
(iso-propyl alcohol)



Material No.: 9059-09

---

Test	Specification
------	---------------

---

---

For Microelectronic Use

Packaging Site: Paris Mfg Ctr & DC

  
Jamie Ethier  
Vice President Global Quality

For questions on this Product Specification please contact Technical Services at 855.282.6867 or +1.610.386.1700  
Avantor Performance Materials, LLC  
100 Mansford Rd, Suite 200, Radnor, PA 19087. U.S.A. Phone 610.386.1700