



Product No. F780
Specifications current as of:
Jun 8 2011

TEST	SPECIFICATION
Assay (HSCH ₂ (CHOH) ₂ CH ₂ SH)	min. 99.5 %
Absorbance of a 0.1M Solution (1-cm path vs water) at:	
280 nm	max. 0.06
260 nm	max. 0.40
Appearance	Passes Test
Melting Point	40 - 43 °C.

For Laboratory, Research or Manufacturing Use

Keep material refrigerated between 2-8°C (36-46°F).

Store Under Nitrogen