

Catalog Number: BDH3116-1L
Catalog Description: Ether, Anhydrous
Grade: Absolute, ACS
Date Current Specification Set: 01-JAN-2014

All batches of this product manufactured after this date are tested to this specification. A new certificate will be published on the web in the unlikely event of specification changes.

Characteristics	Specifications
Peroxide (as H ₂ O ₂)	1 ppm max.
Preservative (BHT) Level	Record only
Residue After Evaporation	0.001% max.
Titration Acid	0.0002 meq/g max.
Water (H ₂ O)	0.03% max.
Assay ((CH ₃ CH ₂) ₂ O)	99.0% min.
Carbonyl (as HCHO)	0.001% max.
Color (APHA)	10 max.
Alcohol (CH ₃ CH ₂ OH)	To pass test (0.05%)
Manufacture Date	Record Actual

Signed on behalf of VWR:



Laura Buholca
Technical Manager

Additional Information (optional):