

Product Specification

Meets ACS Reagent Chemical Requirements,

Test	Specification
Appearance (clear, colorless liquid)	Passes Test
Infrared Spectrum (Conforms to Standard)	Passes Test
ACS - Assay (CH ₃ CHOHCH ₃)	≥ 99.5 %
ACS - Color (APHA)	≤ 10
ACS - Residue after Evaporation	≤ 0.001 %
ACS - Solubility in H ₂ O	Passes Test
ACS - Titrable Acid or Base (meq/g)	≤ 0.0001
ACS - Water (H ₂ O)(by Karl Fischer titrn)	≤ 0.2 %
Carbonyl Compounds - Propionaldehyde	≤ 0.002 %
Carbonyl Compounds - Acetone	≤ 0.002 %

For Laboratory, Research, or Manufacturing Use

Date Received: / /

Date Opened: / /

Packaging Site: Paris Mfg Ctr & DC



Jamie Ethier
Vice President Global Quality