

BOTANICAL **ETHANOL EXTRACTION**

Jacketed reactor base system designed for ethanol extraction and winterization processes. Contact our technical support staff to complete the reactor components, stir motor, chiller, etc. that best suit your needs.

- Designed specifically to operate at very low temperatures
- Configured for an explosion-proof motor
- Configure any of our kilo-scale reactors for winterization
- 10L to 150L

Base Reactor System Includes:

- Mobile support stand
- Jacketed glass vessel
- Head and lid
- Bearing
- Agitators
- Stir shaft
- Clamps
- Bottom outlet valve is equipped with EPDM o-rings

XP Series Stirring

- Hazardous duty motor
- Digital control box
- CE and cULus rated



Also Available...

Contact us for all your reactor needs, for both large and bench scale applications.

KILO-SCALE



Shown with Motor and Control, not included.

BOTANICAL **DISTILLATION KIT**

Add-on distillation kit for use with any Ace kilo scale reactor with our Universal Stand. Kits include a side shelf to attach to your existing stand, receiving flask, distillation splitter, condenser, tubing and adapters. Please contact Ace Technical Support for customization options.

Distillation Add-on Package:

- Side shelf
- 22L spherical receiving flask
- Offset adapter
- Splitter
- Condenser
- Tubing
- Fittings

Chillers and Recirculators

- Lauda
- Polyscience
- Julabo



KILO-SCALE



Shown on 12855 Series Base System Reactor with Motor and Control

Also Available...

Contact us for custom applications and accessories.

avantor[™]
delivered by **VWR**[™]

VWR.COM

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada unless otherwise noted, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2022 Avantor, Inc. All rights reserved.