

# Certificate of Quality



Part Number	MFLX31201-14
Lot Number	S00557TBL1
Part Description	ADTPR 3/4IN TC 1/8IN ID STRL
Expiration Date	30MAR2024
Country of Origin	USA

### Inspection & Testing

- It is the responsibility of the user to conduct tests that are deemed necessary to determine the suitability of this product for any particular use.

### Product Specifications & Compliance

- Compound: Polypropylene with fused TPE gasket
- All components intended for fluid contact are certified by testing and/ or attestations provided by approved vendors to meet the below criteria.
- All products are packaged in an ISO Class 7 Cleanroom.
- Gamma irradiated to a minimum of 25 kGy and a maximum of 40 kGy.
- This product is Animal-Derived Component Free.
- The material meets USP Class VI requirements.
- Date of expiration is set to 5 years from the date of manufacture (DOM).
- It is the responsibility of the user to conduct tests that are deemed necessary to determine the suitability of this product for any particular use.
- This product is not approved or intended for, and should not be used for medical, clinical, surgical or other patient oriented applications.

### Storage

- Storage in original packaging, at ambient conditions and away from direct sunlight, is recommended.

### Certificate Approval

- This certificate was issued and approved by the Avantor Fluid Handling Quality and Regulatory Departments and is valid without signature.