



**Item Number** X1017  
**Item** Xylitol, NF  
**CAS Number** 87-99-0  
**Molecular Formula** C<sub>5</sub>H<sub>12</sub>O<sub>5</sub>  
**Molecular Weight** 152.15

Assay (Anhydrous Basis).....98.5-101.0%

MAXIMUM LIMITS

Water..... 0.5%  
Residue on Ignition..... 0.5%  
Heavy Metals (as Pb)..... 0.001%  
Reducing Sugars.....0.2%  
Other Polyols.....2.0%  
Residual Solvents.....To pass test

***Spectrum Chemicals & Laboratory Products***

---

**Corporate Headquarters:**      **East Coast Facility:**  
14422 S. San Pedro St.      755 Jersey Ave.  
Gardena, CA 90248      New Brunswick, NJ 08901  
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843