

Product Specification

Meets ACS Reagent Chemical Requirements,

Test	Specification
Assay (NaCN) (by Ag titrn)	≥ 95.0 %
Chloride (Cl)	≤ 0.15 %
Phosphate (PO ₄)	≤ 0.02 %
Sulfate (SO ₄)	≤ 0.05 %
Sulfide (S)	≤ 0.005 %
Thiocyanate (SCN)	≤ 0.02 %
Iron (Fe)	≤ 0.005 %
Trace Impurities – Lead (Pb)	≤ 5 ppm

For Laboratory, Research, or Manufacturing Use
Meets Reagent Specifications for testing USP/NF monographs

Packaging Site: Paris Mfg Ctr & DC