



FEEL THE SW
DIFFERENCE™

FOCUSED ON ADVANCED CHEMICAL HAND PROTECTION.




PC-079-095-GR
NEOPRENE POWDER-FREE EXAM GLOVES

Previously known as NeoForm® X8

WHEN THE SITUATION CALLS FOR HIGH PERFORMANCE, OUR NEOPRENE GLOVES PROVIDE SUPERIOR RESISTANCE TO CHEMICALS, PUNCTURES AND ABRASIONS WITH A SECURE GRIP FOR HANDLING SMALL PARTS AND SLIPPERY OBJECTS.

- Pure neoprene prevents latex allergic reactions
- Textured surface delivers excellent grip
- Ergonomic design for comfort
- Optimal thickness for durability and tactile sensitivity
- Powder-free to prevent contamination
- Tear-resistant beaded cuff eases donning
- Tested for use with chemotherapy drugs
- SHA Dermatologically Approved
- 1.5 AQL for wet and dry environments



PC **PC-079-095-GR**

WEIGHT	LENGTH	THICKNESS	COLOR
7.9g (large size)	9.5"	4.7 mil	GREEN



SW Safety Solutions Inc.
33278 Central Avenue, Unit 102, Union City, CA 94587, USA
Tel: +1.510.429.8692 | Fax: +1.510.487.5347
Trademarks and registered trademarks are the property of SW Safety Solutions Inc. and its affiliates.
©2019 SW Safety Solutions Inc. All rights reserved.



TECHNICAL DATA



PC-079-095-GR NEOPRENE POWDER-FREE EXAM GLOVES

Previously known as NeoForm® X8

PRODUCT DESCRIPTION

SPECIFICATIONS

PRODUCT	PowerChem
GRADE	Powder-Free Exam
MATERIAL	Neoprene
THICKNESS	FINGER 6.2mil 0.16mm
	PALM 4.7mil 0.12mm
	CUFF 3.8mil 0.10mm
COLOR	Green
SURFACE	Textured
LENGTH	9.5" 242mm
WEIGHT	7.9g
FIT TYPE	Ambidextrous
CUFF TYPE	Beaded
WATERTIGHT AQL	1.5
ORIGIN COUNTRY	Indonesia

PHYSICAL PROPERTIES

	PRE-AGING	POST-AGING
TENSILE STRENGTH:	18 MPA	18 MPA
ELONGATION:	750%	700%



* Specifications based on size Large glove

STANDARDS & CERTIFICATIONS

EN 420-2003
EN 455-2000-2015
ASTM D6319-10
ASTM F1671-13
ASTM F739-12
REACH EU 1907/2006
EU 2016/425

CAT III COMPLEX DESIGN
FDA 21 CFR 177.2600
Component materials
comply with FDA
regulations for food
contact



Please visit swsafety.com/symbols for information on standards used on this document.

ORDERING

ITEM NO.	SIZE	PALM WIDTH	PACKAGING
K001401	XS 5½ – 6	3.0" 75mm	100 gloves per box 10 boxes per case
K001402	S 6½ – 7	3.4" 86mm	100 gloves per box 10 boxes per case
K001403	M 7½ – 8	3.8" 97mm	100 gloves per box 10 boxes per case
K001404	L 8½ – 9	4.2" 107mm	100 gloves per box 10 boxes per case
K001405	XL 9½ – 10	4.5" 114mm	100 gloves per box 10 boxes per case

APPLICATIONS

- Chemical Handling
- Janitorial Maintenance
- Light Chemical Mixing
- Painting

Keep out of sunlight. Store in a cool, dry place.
Keep away from sources of ozone or ignition.



R_K00140X1119



SW Safety Solutions Inc.
33278 Central Avenue, Unit 102, Union City, CA 94587, USA
Tel: +1.510.429.8692 | Fax: +1.510.487.5347

Trademarks and registered trademarks are the property of SW Safety Solutions Inc. and its affiliates.
©2019 SW Safety Solutions Inc. All rights reserved.

