

Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® Vinyl Nitrile Hybrid Examination Gloves, Blue color, Length 9.5", Powder free, Beaded cuff style
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2022-02-09

VWR Catalog Number	Glove size
76456-892	Small
76456-894	Medium
76456-896	Large
76456-898	X-Large

Quality System Compliance

Products are manufactured under the **ISO 9001:2016 and ISO 13485:2016** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Accuracy:	Watertight test:	Inspection level = G1; AQL 1.5
	Visual Inspection:	Inspection level = G1; AQL 2.5 (major defects) AQL 4.0 (minor defects)
	Dimensional:	Inspection level = S2; AQL 4.0
	Physical Property:	Inspection level = S2; AQL 4.0
	Powder Content:	Inspection level = 5 pieces; Average $\leq 2\text{mg/glove}$
	Glove count:	Inspection level = 5 dispensers; Tolerance 100 ± 2 pieces

Product Specifications

Material	Synthetic Polyvinylchloride (PVC). This product conformed to ASTM D 6319 standards.
Sterilization	Not applicable.
ATP Assay	Not applicable.

DNase & RNase Free	Not applicable.
BSE/TSE	Not applicable.
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	The product is Latex free.
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Please see RoHS information below.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr ⁶⁺), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP), DEHP (di(2- ethylhexyl) phthalate)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2021-07-08 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.
Storage Conditions	Stored in dark and room temperature environment (5-38°C). Protect gloves against ultraviolet light sources, such as sunlight and oxidizing agents. Keep dry.
Shelf Life	5 years.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.