

Product Name:	PRF122WWW/0CAD
Description:	<p>Corepoint Scientific™ White Diamond Series Refrigerator &amp; Freezer Combination with Controlled Auto Defrost Freezer  12 Cu. Ft. capacity (10.5 Ref./1.0 Freezer)  Two exterior swing doors (solid door refrigerator/solid door freezer), right hinged, not reversible  Microprocessor temperature controller, refrigerator and freezer Digital temperature display with sample simulated probe, refrigerator and freezer  Adjustable operating temperature range:  Refrigerator: 36°F - 46°F (2°C - 8°C)  Freezer: 5°F – -13°F (-15°C – -25°C)  Defrost:  Auto defrost cycle every 6 hours: (at Set Point -11°F /-24°C)  Max air temp in empty freezer: 5°F (-15°C)  Max air temp in loaded freezer: 3°F (-16°C)  Max sample temp in empty freezer: 0°F (-18°C)  Max sample temp in loaded freezer: -4°F (-20°C)  Refrigerator/Freezer freestanding thermometer with 3 year certificate of calibration, with min/max memory, °F/°C switchable Audible and visual high and low temperature alarms, refrigerator and freezer  Remote alarm contacts, refrigerator and freezer  Six adjustable shelves with guard rail on back to prevent product from being pushed against back wall, Refrigerator  Two adjustable shelves with guard rail on back to prevent product from being pushed against back wall, Freezer  Bottom footed shelf to allow air circulation, Refrigerator and Freezer  Vaccine storage power cord warning label  Forced draft circulation Refrigerator  Keyed door locks  Magnetic door gasket for positive seal  Integrated door handles  Leveling legs  Stacking kit  LED interior light (Refrigerator)  Two probe access ports (3/8"); one, refrigerator/one, freezer  White powder coated exterior  HFC-free refrigerant (R600a), foam insulation and packaging  Pharmacy refrigerator/freezer toolkit and temperature logs  Three year parts and labor warranty, plus an additional two year compressor parts warranty  Pyxis®, Omnicell® and AcuDose RX® compatible  Exterior dimensions: 23 ¾"W x 26 ½" D x 80 ¼"H  Controlled Auto Defrost freezer, Patented Technology  Cycle defrost refrigerator</p>

	<p>115V, 60 Hz, 2.3 Amps refrigerator, 1 Amp freezer, 1/3 HP refrigerator, 1/8 HP freezer UL/C-UL listed Shipping weight: 316 lbs. <b>In stock, with next day shipping</b> *Please note this unit requires 2 outlets.</p>
--	---