

# Sterile TechniSat®

TechniCloth® wipers pre-wetted with 70% isopropyl alcohol and sterilized



## Description

Sterile TechniSat® is constructed from a hydroentangled polyester/cellulose blend that is pre-wetted with 70% isopropyl alcohol and packaged in a resealable flexpack. Each package is gamma irradiated to 10<sup>-6</sup> Sterility Assurance Level.

Sterile TechniSat® is used for alcohol wiping in sterile environments in pharmaceutical aseptic fill areas, biotechnology manufacturing facilities, microbiological laboratories, sterile suites and prep rooms.

## Features

- Gamma irradiated to 10<sup>-6</sup> Sterility Assurance Level, according to AAMI Guidelines ensures sterility
- Lot-specific information simplifies record keeping:
  - Certificate of Irradiation confirms radiation dosage,
  - Certification of Compliance, and
  - Expiration date marked on each package
- Pre-wetted with 70% isopropyl alcohol/30% deionized water, filtered through a 0.2µm filter
- Binder-free, hydroentangled polyester/cellulose substrate
- Independent quality control audits for sterility assurance

## Benefits

- Optimizes solution usage and reduces spill potential
- Eliminates problems associated with filtration, sterilization and mixing of alcohol and DI water
- Packaging permits alcohol wipedown of exterior bag before introduction into sterile suites
- Offers low levels of solvent extractables and particle generation

## Applications

- Designed for cleaning and polishing equipment and environmental surfaces during and following the production flow
- Ideal for wiping gloved hands
- Suitable for sterile wiping of production equipment and laminar flow hoods requiring low solvent extractables and low particle generation

## Products

Number	Description	Packaging
<b>TX3214</b>	Sterile TechniSat® 9" x 11" (23 cm x 28 cm) pre-wetted TechniCloth® wipers	50 wipers/flexpack 24 flexpacks/case
<b>TX3217</b>	Sterile TechniSat® 9" x 11" (23 cm x 23 cm) per-wetted TechniCloth® wipers	20 wipers/flexpack 24 flexpacks/case

Custom cleaning solutions available upon request.

**North America**  
1210 South Park Drive  
Kemersville, NC 27284  
Tel 336.996.7046  
[www.texwipe.com](http://www.texwipe.com)  
[info@texwipe.com](mailto:info@texwipe.com)

**Europe/Middle East**  
Saffierlaan 5  
2132 VZ Hoofddorp  
The Netherlands  
Tel +31 88 1307 400

**Asia/Pacific**  
50 Tagore Lane  
#02-01 Entrepreneur Centre  
Singapore 787494  
Tel +65 6468 9433