

Faster. Safer. More efficient harvesting. Now it's all in the bag.



1 L HarvestLine Bottle Liners

Safe, efficient workflows are as critical to your study as quality samples. Now you can harvest high-quality samples in fewer steps from prep to cleanup. With new HarvestLine sterile bottle liners, optimized speed, safety and pelleting efficiency are in the bag.

Optimize lab safety and efficiency

- Prevent contamination
- Reduce exposure to potentially hazardous biologics
- Reduce repetitive motions
- Prevent liquid spillage
- Reduce needed storage space

Accelerate workflows

- Achieve maximum sample yield
- Minimize non-value-added tasks e.g., bottle cleaning, autoclaving, scraping
- Expedite freeze/thaw lysing



It's all in the **bag**

1 L HarvestLine System Liners



Ordering Information

VWR Cat. No.	Description	
Liners		
76292-060	1 L HarvestLine System Liner, qty 6 (non-sterile)	Material: Polyetherurethane, USP 23 Class VI
76292-062	1 L HarvestLine System Liner, qty 6 (sterile w/peel pack)	
BK369256	1 L HarvestLine System Liner, qty 120 (non-sterile)	
Bottles		
76294-186	1 L Polypropylene Wide-mouth Bottle Assembly (set of 2)	
76294-184	1 L Polycarbonate Wide-mouth Bottle Assembly (set of 2)	
Rotors		
BK363688	JLA-8.1000 Rotor – 6 x 1 L; 8,000 rpm; 15,970 x g	
BK366754	JLA-9.1000 Rotor – 4 x 1 L; 9,000 rpm; 16,800 x g	
Accessories		
76292-132	1 L HarvestLiner Kit (includes rack and funnel)	

Visit vwr.com/beckman or contact your VWR Sales Representative for more information.

800 932 5000 | VWR.COM

Prices and product details are current when published and subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2019 Avantor, Inc. All rights reserved.

