

Sodium Hydroxide, 10N Volumetric Solution
BAKER ANALYZED® Reagent



Material No.: 5674-02

Revision No.: 0

Product Specification

Test	Specification
Appearance (Clear, colorless liquid)	Passes Test
Normality	9.9 - 10.1
Sodium Carbonate (Na_2CO_3)	$\leq 2.0 \%$
Chloride (Cl)	$\leq 5 \text{ ppm}$
ACS - Heavy Metals (as Pb)	$\leq 1 \text{ ppm}$
Trace Impurities - Iron (Fe)	$\leq 1 \text{ ppm}$

For Laboratory, Research, or Manufacturing Use

Standardization at 25°C traceable to NIST Standard Reference Material.

SRM No

IMPORTANT: Material will freeze if stored below 41°F (5°C).

Storage Condition: IMPORTANT: Material will freeze if stored below 41 °F (5°C). Protect from air to avoid absorption of carbon dioxide.

Packaging Site: Phillipsburg Mfg Ctr & DC


Jamie Ethier
Vice President Global Quality

For questions on this Product Specification please contact Technical Services at 855.282.6867 or +1.610.386.1700

Avantor Performance Materials, LLC

100 Matsonford Rd, Suite 200, Radnor, PA 19087. U.S.A. Phone 610.386.1700