

Sodium Sulfate, Anhydrous  
U.S.P.



Material No.: 8028-20

Revision No.: 0

## Product Specification

Meets U.S.P Requirements, GMP Manufactured Product

Test	Specification
USP - Assay (Na <sub>2</sub> SO <sub>4</sub> ) (dried basis)	≥ 99.0 %
USP - Acidity or Alkalinity (mL)	≤ 0.50
USP - Chloride (Cl)	≤ 0.005 %
USP - Identification A	Passes Test
USP - Identification B	Passes Test
USP - Insoluble Matter	≤ 0.01 %
USP - Loss on Drying at 105°C	≤ 0.5 %
Thru U.S. No. 50 Sieve, For Information Only	%

Bulk Pharmaceutical Chemical

GMP Manufactured Product

CAUTION: For Manufacturing, processing or repackaging

No Class 1,2,3 or other solvents are used or produced in the manufacturing or purification of the product.

Metallic Residues: No metal catalysts or metal reagents, as defined by EMA Guideline EMEA/CHMP/SWP/4446/2000, are used in the production of this material.

Appearance - Colorless to slightly colored crystalline solid.

Packaging Site: Paris Mfg Ctr & DC

  
Jamie Ethier  
Vice President Global Quality

For questions on this Product Specification please contact Technical Services at 855.282.6867 or +1.610.386.1700

Avantor Performance Materials, LLC

100 Matsonford Rd, Suite 200, Radnor, PA 19087. U.S.A. Phone 610.386.1700