



Toluene
BAKER ANALYZED ACS
Reagent

Product No. 9460
 Specifications current as of: Jun 8
 2011

TEST	SPECIFICATION
Meets A.C.S. Specifications	
Meets Reagent Specifications for testing USP/NF monographs	
Assay (C ₆ H ₅ CH ₃) (by GC)	min. 99.5 %
Color (APHA)	max. 10
Residue after Evaporation	max. 0.001 %
Substances Darkened by H ₂ SO ₄	Passes Test
Sulfur Compounds (as S)	max. 0.003 %
Water (by KF, coulometric)	max. 0.03 %

Product Information (not specifications):	
Boiling Point (typical)	°C. 110.6

For Laboratory, Research or Manufacturing Use