



Product Specification

Catalog Number: BDH9280

Catalog Description: BDH SODIUM BICARB ACS

Grade: ACS

Date Current Specification Set: NOV. 2012

All batches of this product manufactured after this date are tested to this specification. A new certificate will be published on the web in the unlikely event of specification changes.

Characteristics	Specifications
Ammonium	$\leq 0.0005\%$
Calcium	$\leq 0.02\%$
Chloride	$\leq 0.003\%$
Heavy Metals (as Pb)	$\leq 0.0005\%$
Insolubles	$\leq 0.015\%$
Iron	$\leq 0.001\%$
Magnesium	$\leq 0.005\%$
Phosphate	$\leq 0.001\%$
Potassium	$\leq 0.005\%$
Purity	99.7 - 100.3 %
Sulfur Compounds	$\leq 0.003\%$

Signed on behalf of VWR:

Robin Crocker
Quality Control Manager