

glove liners and inspection gloves

HIGHLY REUSABLE NYLON LINER

- 100% continuous nylon filament makes the liner thin and comfortable
- Fully launderable up to 30 times without shrinkage, reducing cost of waste disposal
- · Available in full and half finger
- Ambidextrous





Size	Part No.	Part No.
W	M006W	32935-566
M	M006M	32935-568
*Sold by the dozen pair		

WLI

CUT RESISTANT LINER

- · Stainless steel and polyester fiber construction provides cut resistance
- Highly launderable
- Use as a liner or stand alone
- · Black color hides dirt and grime
- Ambidextrous
- ANSI Level 3
- · Sold by the dozen pair

Size	WLI Part No.	VWR Part No.
S	M281S	32935-392
M	M281M	32935-394
L	M281L	32935-396
XL	M281XL	32935-398



KEVLAR® LINER

- Lightweight
- · Highly economical
- Offers minor cut resistance
- Ambidextrous
- ANSI Level 2
- Sold by the dozen pair

Size	WLI Part No.	VWR Part No.
S	M79S	80091-677
M	M79M	32935-155
L	M79L	32935-156



REUSABLE KNITTED NYLON LINER

- Stretch nylon construction
- Form fitting
- Excellent when a high level of dexterity is required
- Ambidextrous
- · Sold by the dozen pair

Full Fin	ger		Half Fin	ger		
Size	WLI Part No.	VWR Part No.	Size	WLI Part No.	VWR Part No.	
OS	M089	32935-535	OS	M088	32935-479	

CUT RESISTANT LINER

- High performance fiber construction
- · High reusability
- Stretchable
- Ambidextrous
- ANSI Level 2
- · Sold by the dozen pair

Size	WLI Part No.	VWR Part No.
XS	M214XS	32935-509
S	M214S	32935-481
M	M214M	32935-483
L	M214L	32935-485
XL	M214XL	32935-487



INSPECTION GLOVE

- Light 100% cotton liner
- Single use
- · Superior tactile sensitivity
- · Available in light and medium





е	WLI Part No.	VWR Part No.	Size	WLI Part No.	VWR Part No.
	Y6701W Y6701L	32935-461 32935-460	W	Y6703W Y6703L	32935-550 32935-552

medical liners

SPECTEC® STERILE LINER

- 100% high performance fiber construction provides cut resistance
- Excellent tactile sensitivity
- Sterile
- Ambidextrous
- Sold by the box, 20 pieces/box
- ANSI Level 1

Size	WLI Part No.	VWR Part No.
S	M102S	80091-126
M	M102M	80091-124
L	M102L	80091-122



NYLON LINER

- 100% continuous nylon filament
- Used to protect hands from prolonged exposure to latex
- Fully launderable up to 20 times without shrinking
- Ambidextrous
- Non-sterile
- Sold by the box, 25 pieces/box

Extende	ed Cuff, Full Fin	ger	Regula	r Cuff, Full Fing	er
Size	WLI Part No.	VWR Part No.	Size	WLI Part No.	VWR Part No.
S M	M113S M113M	32935-100 32935-102	S M	M115S M115M	14224-270 13500-888
1	M113I	32035-104	1	M115I	1/100/1-070

^{* 20} pairs per Poly bag

leather

GRAIN COWHIDE

- Premium quality cowhide leather
- Gunn cut
- Straight thumb
- · Double shirred elastic wrist
- Bound hem
- Sold by the dozen pair

Size	WLI Part No.	VWR Part No.
S	Y01225S	32935-448
M	Y01225M	32935-449
L	Y01225L	32935-450
XL	Y01225XL	32935-451
XXL	Y01225XXL	89030-818



STANDARD SHOULDER SPLIT **LEATHER PALM**

- Standard shoulder split palm
- Wing thumb for increased grip, comfort and flexibility
- Safety cuff
- Economy choice
- · Sold by the dozen pair

Size	WLI Part No.	VWR Part No.
S	Y3407S	80092-126
L	Y3407L	80092-124



GRAIN GOATSKIN

- · Quality goatskin leather
- Excellent wear and abrasion
- · Gunn cut for increased wear life
- Keystone thumb for extra comfort
- Bound hem
- Sold by the dozen pair

Size	WLI Part No.	VWR Part No.
S	Y0107S	80091-856
M	Y0107M	80091-854
L	Y0107L	80091-852
XL	Y0107XL	80091-850



SELECT SHOULDER SPLIT LEATHER PALM

- Selected shoulder palm
- Full leather palm
- Knuckle strap
- Wing thumb
- Sold by the dozen pair



VWR

Y3107I

Y3107XL

Part No.

80092-092 80092-090

80092-088

80092-086

Safety Cuff	
-------------	--

Y3106L

Y3106XL

Size

S

М

XL

WLI Part No.	VWR Part No.	Size	WLI Part No.	
Y3106S Y3106M	80092-082	S	Y3107S	

ıff	Gaunt

80092-078

80092-076

heat

WOOL-LINED KEVLAR®

- Protects up to 700°F
- 100% Kevlar outer shell, which provides excellent thermal properties
- Fully lined with virgin wool
- Glove length 13"
- · Loop-out, clute pattern
- Sold by the dozen pair

Size	WLI Part No.	VWR Part No.
ΧI	305KWI	32935-691



KEVLAR® DUCK GAUNTLET

- Protects up to 700°F
- 100% Kevlar® outer shell insulates against heat
- Fully lined with virgin wool
- 5" flame resistant duck gauntlet cuff
- · Loop-out, clute pattern
- · Sold by the dozen pair

Size	WLI Part No.	VWR Part No.
XL	625	80091-235



AUTOCLAVE GLOVE

- Protects to 250°F
- · Loop-out terrycloth
- Clute cut
- · Gauntlet cuff
- Sold by the dozen pair

Cuff Size	WLI Part No.	VWR Part No.
5"	422-5	80091-460
11"	422-11	80091-458



BROWN & WHITE TERRY GLOVE

- Heat protection up to 450°F
- Extra heavy weight terry cloth
- Loop-out
- Gunn Cut
- · Sold by dozen pair

Cuff	WLI	VWR
Size	Part No.	Part No.
Safety	644HRL	80091-490
Gauntlet	636HRL	32935-698



palm dipped

KEVLAR® WITH FOAM NITRILE PALM

- Light weight stretch Kevlar with black foam nitrile palm coating
- Stretchy Kevlar/Lycra® blend provides high level of dexterity and comfort
- Foam nitrile coating provides grip and resistance to abrasions
- ANSI Level 2
- · Sold by the dozen pair

Size	WLI Part No.	VWR Part No.
XS	Y9256XS	89030-808
S	Y9256S	89092-390
M	Y9256M	89092-388
L	Y9256L	89092-386
XL	Y9256XL	89092-392



HDPE WITH POLYURETHANE PALM

- High performance polyethylene cut resistant fiber (HPPE)
- Polyurethane palm coating
- Exceptional comfort, dexterity, cut resistance and grip
- ANSI Level 2
- Sold by the dozen pair

Size	WLI Part No.	VWR Part No.
XS S M	Y9265XS Y9265S Y9265M	89030-788 89092-790 89092-792
L XL	Y9265L Y9265XL	89092-792 89092-794 89092-796



NYLON WITH POLYURETHANE PALM

- Nylon with polyurethane palm
- Coated palm and fingertips provide additional protection and grip
- Seamless design offers increased dexterity and flexibility
- Sold by the dozen pair

Size	WLI Part No.	VWR Part No.
XS	Y9277XS	89030-798
S	Y9277S	89092-800
M	Y9277M	89092-802
L	Y9277L	89092-804
XL	Y9277XL	89092-806



POLY/COTTON KNIT WITH LATEX PALM

- Polyester/cotton string knit shell
- Latex coating offers superior grip and abrasion resistance
- · Flexible and comfortable
- · Gray color conceals dirt and grime
- · Sold by the dozen pair

Size	WLI Part No.	VWR Part No.
XS	Y9240XS	89030-810
S	Y9240S	89092-352
M	Y9240M	89092-350
L	Y9240L	89092-348
XL	Y9240XL	89092-346



cut

SILVER TALON®

- Unmatched combination of cut resistance, dexterity and durability
- Second skin fit allows for increased flexibility and comfort
- Silver color better conceals stains
- Ambidextrous
- ANSI Level 5
- Sold by the piece

Size	WLI Part No.	VWR Part No.
XXS	134524	14231-960
XS	134525	32935-158
S	134526	32935-160
M	134527	32935-162
L	134528	32935-164
XL	134529	32935-166



HANDGUARD® II

- Heavy duty patented combination of high performance cut resistant fibers and stainless steel
- Unmatched combination of comfort and cut resistance
- Ambidextrous
- ANSI Level 5
- · Sold by the piece

Size	WLI Part No.	VWR Part No.
XS	333019	32935-178
S	333021	32935-138
M	333023	32935-140
L	333025	32935-142
XL	333027	32935-182





Prices and product details are current when published; subject to change without notice. [Certain products may be limited by federal, state, provincial, or local regulations, IWR makes no claims or warranties concerning sustainabled/green products. Any claims concerning sustainabled/green products are the sole claims of the manufacturer and not those of VMR International, LLC, All prices are in US dollars unless otherwise noted, Offers valid in US, void where prohibited by law or company policy, while supplies last, IVMR, the VMR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks, of VMR International, LLC and Its related companies. All other marks referenced are registered by their respective owner(s), I Visit vwx.com to view our privacy policy, trademarks owners and additional disclaimers, © 2012 VWR International, LLC, All rights reserved.