

Name of the Product	: Uridine Cell Culture Tested
Code No.	: TC181
Section 1	: Chemical Identification Code No. : TC181 Name of the Product : Uridine Cell Culture Tested Produced by : HiMedia Laboratories Pvt. Ltd. Address : 23, Vadhani Indl. Estate, LBS Marg, Mumbai 400 086, India. Tel. No. : 2500 0970, 2500 1607 Fax No. 022 2500 2468
Section 2	: Uridine CAS No : 58-96-8
Section 3	: Hazards Identification Hazard : Not classified as hazardous.
Section 4	: First - Aid Measures No specific measures necessary.
Section 5	: Fire Fighting Measures Not combustible.
Section 6	: Accidental Release Measures No specific measures necessary.
Section 7	: Handling and Storage Handling - Refer to Section 8 Storage - Store at 2 - 8°C
Section 8	: Exposure Controls/Personal Protection Wear appropriate NIOSH/MSHA-approved respirator, chemical-resistant gloves, safety goggles, other protective clothing. Mechanical exhaust required.
Section 9	: Physical and Chemical Properties Appearance : Amorphous powder Colour : White
Section 10	: Stability and Reactivity Stability : Product is stable if stored as per the conditions specified under storage of Section No. 7. Product loses its potency/performance above 45°C. Conditions to avoid : Heat and light.

Section 11 : Toxicological Information

Non toxic.

Section 12 : Ecological Information

Data not available

Section 13 : Disposal Considerations

No special disposal method required except that it be in accordance with current and local authority regulation.

Section 14 : Transport Information

UN No. : Not applicable.

Section 15 : Regulatory Information

Risk Phrases : Not applicable

Safety Phrases : Avoid contact with skin and eyes

Section 16 : Other Information

The information contained in this data sheet represents the best information currently available to us. However, no warranty is made with respect to its completeness and we assume no liability resulting from its use. The information is offered solely for user's obligation to investigate and determine the suitability of the information for their particular purpose.

