



### VWR-SLD BSKT-FULL-5.2\*



Compatible Refrigerator Size	5.2 ft <sup>3</sup>
Max Quantity	3
Basket type	Sliding
Width x Depth x Height	16 7/8" x 13 5/8" x 3"
Max Load	30 lbs.


\*Factory install only

### VWR-BSKT-BASE-FIXED-5.2



Compatible Refrigerator Size	5.2 ft <sup>3</sup>
Max Quantity	1
Basket type	Fixed
Width x Depth x Height	18 3/8" x 8 1/2" x 6 1/4"
Max Load	44 lbs.


### VWR-SLD BSKT-FULL-10.5\*



Compatible Refrigerator Size	10.5 ft <sup>3</sup>
Max Quantity	6
Basket type	Sliding
Width x Depth x Height	16 7/8" x 15 3/4" x 3"
Max Load	35 lbs.


\*Factory install only

### VWR-BSKT-BASE-FIXED-10.5



Compatible Refrigerator Size	10.5 ft <sup>3</sup>
Max Quantity	1
Basket type	Fixed
Width x Depth x Height	18 3/8" x 10 3/8" x 6 1/4"
Max Load	55 lbs.


### VWR-SLD BSKT-FULL-5.2



HCUCFS0504 / 76300-508  
 HCUCFS0504G / 10819-876  
 †VWR-HC-RFC7 / 76580-180  
 †HCRFC7S / 76488-190  
 †HCRFC7SA / 76488-192  
 †HCRFC1020 / 76488-178

†HCRFC1020G / 76488-184  
 †HCRFC1030 / 76488-180  
 †HCRFC1030G / 76488-186  
 †HCRFC1040 / 76488-182  
 †HCRFC1040G / 76488-188

### VWR-BSKT-BASE-FIXED-5.2



HCUCBI0404 / 10819-880\*  
 HCUCBI0404G / 10819-882\*  
 HCUCBI0404SS / 10819-886\*  
 HCUCBI0404SAD / 76422-930\*  
 HCUCBI0404AD / 10819-884\*  
 HCUCFS0504 / 76300-508  
 HCUCFS0504G / 10819-876  
 †VWR-HC-RFC7 / 76580-180

†HCRFC7S / 76488-190  
 †HCRFC7SA / 76488-192  
 †HCRFC1020 / 76488-178  
 †HCRFC1020G / 76488-184  
 †HCRFC1030 / 76488-180  
 †HCRFC1030G / 76488-186  
 †HCRFC1040 / 76488-182  
 †HCRFC1040G / 76488-188


### VWR-SLD BSKT-FULL-10.5



VWR-HC-10PG / 76580-148  
 VWR-HC-10PS / 76580-150

†VWR-HC-RFC12G / 76580-178  
 †VWR-HC-RFC12 / 76580-176

### VWR-BSKT-BASE-FIXED-10.5



VWR-HC-10PG / 76580-148  
 VWR-HC-10PS / 76580-150

†VWR-HC-RFC12G / 76580-178  
 †VWR-HC-RFC12 / 76580-176

\*Model available in left hinge

†Compatible with refrigerator compartment only