



Item Number F1135
Item Formic Acid, FCC
CAS Number 64-18-6
Molecular Formula CH₂O₂
Molecular Weight 46.03

Assay (CH₂O₂)..... Min. 85.0%
Dilution Test..... To pass test

MAXIMUM LIMITS

Acetic Acid.....0.4%
Sulfate (SO₄).....0.004%

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: **East Coast Facility:**
14422 S. San Pedro St. 755 Jersey Ave.
Gardena, CA 90248 New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843