



MATERIAL SAFETY DATA SHEET

according to the Global Harmonized System

Revision Date 02/04/2013

Version 1.1

SECTION 1. Identification

Product identifier

Product number	102083
Product name	Calcium chloride fused, granular about 0.5-2.0 mm for elementary analysis

Relevant identified uses of the substance or mixture and uses advised against

Identified uses	Reagent for analysis
-----------------	----------------------

Details of the supplier of the safety data sheet

Company	EMD Millipore Corporation 290 Concord Road, Billerica, MA 01821, United States of America SDS Phone Support: +1-978-715-1335 General Inquiries: +1-978-751-4321 Monday to Friday, 9:00 AM to 4:00 PM Eastern Time (GMT-5) e-mail: mm_sds@merckgroup.com
---------	--

Emergency telephone	613-996-6666 CANUTEC (Canada) +1-703-527-3887 CHEMTREC (International) 24 Hours/day; 7 Days/week
---------------------	--

SECTION 2. Hazards identification

GHS Classification

Eye irritation, Category 2, H319

For the full text of the H-Statements mentioned in this Section, see Section 16.

GHS-Labeling

Hazard pictograms



Signal Word
Warning

Hazard Statements
H319 Causes serious eye irritation.

Precautionary Statements

MATERIAL SAFETY DATA SHEET
according to the Global Harmonized System

Product number 102083 Version 1.1
Product name Calcium chloride fused, granular about 0.5-2.0 mm for elementary analysis

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Other hazards

None known.

SECTION 3. Composition/information on ingredients

Formula	CaCl ₂	CaCl ₂ (Hill)
CAS-No.	10043-52-4	
Molar mass	110.98 g/mol	

Hazardous ingredients

Chemical Name (Concentration)

CAS-No.

calcium chloride (>= 90 % - <= 100 %)

10043-52-4

SECTION 4. First aid measures

Description of first-aid measures

Inhalation

After inhalation: fresh air.

Skin contact

After skin contact: wash off with plenty of water. Remove contaminated clothing.

Eye contact

After eye contact: rinse out with plenty of water. Immediately call in ophthalmologist.

Ingestion

After swallowing: immediately make victim drink water (two glasses at most). Consult a physician.

Never give anything by mouth to an unconscious person.

Most important symptoms and effects, both acute and delayed

irritant effects, Stomach/intestinal disorders

Indication of any immediate medical attention and special treatment needed

No information available.

SECTION 5. Fire-fighting measures

Extinguishing media

Suitable extinguishing media

Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

Unsuitable extinguishing media

For this substance/mixture no limitations of extinguishing agents are given.

Special hazards arising from the substance or mixture

MATERIAL SAFETY DATA SHEET
according to the Global Harmonized System

Product number 102083 Version 1.1
Product name Calcium chloride fused, granular about 0.5-2.0 mm for elementary analysis

Not combustible.
Ambient fire may liberate hazardous vapors.
Fire may cause evolution of:
Hydrogen chloride gas

Advice for firefighters

Special protective equipment for fire-fighters

Stay in danger area only with self-contained breathing apparatus. Prevent skin contact by keeping a safe distance or by wearing suitable protective clothing.

Further information

Suppress (knock down) gases/vapors/mists with a water spray jet. Prevent fire extinguishing water from contaminating surface water or the ground water system.

SECTION 6. Accidental release measures

Personal precautions, protective equipment and emergency procedures

Advice for non-emergency personnel: Avoid inhalation of dusts. Avoid substance contact. Ensure adequate ventilation. Evacuate the danger area, observe emergency procedures, consult an expert.

Advice for emergency responders: Protective equipment see section 8.

Environmental precautions

Do not empty into drains.

Methods and materials for containment and cleaning up

Cover drains. Collect, bind, and pump off spills.

Observe possible material restrictions (see sections 7 and 10).

Take up dry. Dispose of properly. Clean up affected area. Avoid generation of dusts.

SECTION 7. Handling and storage

Precautions for safe handling

Observe label precautions.

Conditions for safe storage, including any incompatibilities

Tightly closed. Dry.

Storage temperature: no restrictions.

SECTION 8. Exposure controls/personal protection

Exposure limit(s)

Ingredients

Basis	Value	Threshold limits	Remarks
-------	-------	------------------	---------

calcium chloride 10043-52-4

CAD ON OEL	Time Weighted Average (TWA _{EV}):	5 mg/m ³	
------------	---	---------------------	--

Engineering measures

MATERIAL SAFETY DATA SHEET
according to the Global Harmonized System

Product number 102083 Version 1.1
Product name Calcium chloride fused, granular about 0.5-2.0 mm for elementary analysis

Technical measures and appropriate working operations should be given priority over the use of personal protective equipment.

Individual protection measures

Protective clothing should be selected specifically for the workplace, depending on concentration and quantity of the hazardous substances handled. The chemical resistance of the protective equipment should be inquired at the respective supplier.

Hygiene measures

Immediately change contaminated clothing. Apply skin- protective barrier cream. Wash hands and face after working with substance.

Eye/face protection

Safety glasses

Hand protection

Chemical-resistant, impervious gloves complying with an approved standard should be worn at all times when handling chemical products if a risk assessment indicates this is necessary.

Respiratory protection

required when dusts are generated.

Use a properly fitted, air-purifying or air-fed respirator complying with an approved standard if a risk assessment indicates this is necessary. Respirator selection must be based on known or anticipated exposure levels, the hazards of the product and the safe working limits of the selected respirator.

SECTION 9. Physical and chemical properties

Physical state	solid
Color	white
Odor	odorless
Odor Threshold	not applicable
pH	ca. 8 - 10 at 100 g/l 68 °F (20 °C)
Melting point	772 °C
Boiling point/boiling range	> 2,912 °F (> 1,600 °C)
Flash point	does not flash
Evaporation rate	No information available.
Flammability (solid, gas)	No information available.
Lower explosion limit	not applicable
Upper explosion limit	not applicable
Vapor pressure	No information available.

MATERIAL SAFETY DATA SHEET
according to the Global Harmonized System

Product number 102083 Version 1.1
Product name Calcium chloride fused, granular about 0.5-2.0 mm for elementary analysis

Relative vapor density	No information available.
Relative density	2.15 g/cm ³ at 68 °F (20 °C)
Water solubility	740 g/l at 68 °F (20 °C)
Partition coefficient: n-octanol/water	No information available.
Autoignition temperature	No information available.
Decomposition temperature	No information available.
Viscosity, dynamic	No information available.
Explosive properties	Not classified as explosive.
Ignition temperature	not combustible

SECTION 10. Stability and reactivity

Reactivity

exothermic dissolution process with water

Chemical stability

The product is chemically stable under standard ambient conditions (room temperature) .

Possibility of hazardous reactions

Exothermic reaction with:

boron trifluoride, vinylmethyl ether

Generates dangerous gases or fumes in contact with:

Metals, Zinc, (generation of hydrogen)

Conditions to avoid

Exposure to moisture.

Incompatible materials

no information available

Hazardous decomposition products

in the event of fire: See section 5.

SECTION 11. Toxicological information

Information on toxicological effects

Likely route of exposure

Eye contact, Skin contact, Ingestion

MATERIAL SAFETY DATA SHEET
according to the Global Harmonized System

Product number 102083 Version 1.1
Product name Calcium chloride fused, granular about 0.5-2.0 mm for elementary analysis

Acute oral toxicity

LD50 rat: 1,000 mg/kg (IUCLID)

Acute inhalation toxicity

Symptoms: Possible damages:., mucosal irritations

Acute dermal toxicity

LD50 rat: 2,630 mg/kg
(IUCLID)

Skin irritation

rabbit

Result: No irritation

OECD Test Guideline 404

slight irritation

Eye irritation

rabbit

Result: Eye irritation

OECD Test Guideline 405

Causes serious eye irritation.

Genotoxicity in vitro

Ames test

Result: negative

(Lit.)

Specific target organ systemic toxicity - single exposure

The substance or mixture is not classified as specific target organ toxicant, single exposure.

Specific target organ systemic toxicity - repeated exposure

The substance or mixture is not classified as specific target organ toxicant, repeated exposure.

Aspiration hazard

Regarding the available data the classification criteria are not fulfilled.

Carcinogenicity

IARC	No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.
OSHA	No ingredient of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by OSHA.
NTP	No ingredient of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.
ACGIH	No ingredient of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.

MATERIAL SAFETY DATA SHEET
according to the Global Harmonized System

Product number	102083	Version 1.1
Product name	Calcium chloride fused, granular about 0.5-2.0 mm for elementary analysis	

Further information

Systemic effects:

After swallowing of large amounts:

Stomach/intestinal disorders

Further data:

Handle in accordance with good industrial hygiene and safety practice.

SECTION 12. Ecological information

Ecotoxicity

Toxicity to fish

LC50 *Lepomis macrochirus* (Bluegill sunfish): 10,650 mg/l; 96 h (IUCLID)

Toxicity to daphnia and other aquatic invertebrates

EC50 *Daphnia magna* (Water flea): 144 mg/l; 48 h (IUCLID)

Toxicity to algae

IC50 algae: 3,130 mg/l; 120 h (IUCLID)

Persistence and degradability

No information available.

Bioaccumulative potential

No information available.

Mobility in soil

No information available.

Other adverse effects

Additional ecological information

Discharge into the environment must be avoided.

SECTION 13. Disposal considerations

The information presented only applies to the material as supplied. The identification based on characteristic(s) or listing may not apply if the material has been used or otherwise contaminated. It is the responsibility of the waste generator to determine the toxicity and physical properties of the material generated to determine the proper waste identification and disposal methods in compliance with applicable regulations. Disposal should be in accordance with applicable regional, national and local laws and regulations.

SECTION 14. Transport information

Land transport (DOT)

Not classified as dangerous in the meaning of transport regulations.

Air transport (IATA)

Not classified as dangerous in the meaning of transport regulations.

Sea transport (IMDG)

Not classified as dangerous in the meaning of transport regulations.

MATERIAL SAFETY DATA SHEET
according to the Global Harmonized System

Product number 102083 Version 1.1
Product name Calcium chloride fused, granular about 0.5-2.0 mm for elementary analysis

SECTION 15. Regulatory information

United States of America

Canada

WHMIS Classification

D2B Toxic Material Causing Other Toxic Effects

Eye irritant

This product has been classified according to the hazard criteria of the CPR and the MSDS contains all of the information required by the CPR.

Notification status

TSCA: On TSCA Inventory

DSL: All components of this product are on the Canadian DSL.

SECTION 16. Other information

Training advice

Provide adequate information, instruction and training for operators.

Full text of H-Statements referred to under sections 2 and 3.

H319 Causes serious eye irritation.

Key or legend to abbreviations and acronyms used in the safety data sheet

Used abbreviations and acronyms can be looked up at www.wikipedia.org.

Revision Date 02/04/2013

The information contained herein is based on the present state of our knowledge. It characterizes the product with regard to appropriate safety precautions. It does not represent a warranty of any product properties and we assume no liability for any loss or injury which may result from the use of this information. Users should conduct their own investigations to determine the suitability of the information.

All rights reserved. Millipore and the "M" Mark are registered trademarks of Merck KGaA, Darmstadt, Germany.