

Access DNA/RNA & Proteins from Any Sample Type



Lyse sample in 40 seconds or less!


High Yields

Flexible Sample Holders.

- Ambient to Cryogenic
- 2 mL to 50 mL capacity

FastPrep-24™ 5G

No Sample is too tough for the FastPrep-24 5G!



MP Biomedicals - It's What's Inside That Counts!

NEW FEATURES	BENEFITS
Increased Power Settings to 10 M/S	Highest Power Available Release analytes of interest into protective buffer quicker. Lyse Tougher Samples
Time: 1-120 Seconds Cycles: 1-9 Cycle Pause: 1-300 seconds	Less User Intervention True Walk-Away Operation
Touchscreen User Interface	Intuitive and Easy to Operate
>70 Recommended Programs	Validated, for Quicker Optimizations,
Manual Program-N-Store (Save up to 12 custom programs)	Saves Time and Knowledge
Export Data Function	Archive Assays and Results
USB Port	Easy Software Upgrades
Sound Interface	Audible Program Run Status Alerts
Strobe Effect	Visual Optimization Tool

Lyse Efficiently:

- Tissue, Tumors, Bones
- Bacteria gram + or -
- Soil
- Yeast, Fungi, Spores
- Plants, Seeds, Roots
- Feces

Cat. No.	Description
75784-596	FastPrep-24 5G Tissue Homogenizer

Optimize your research with
MP Lysing Matrix and FastDNA® Purification Kits
for the most intact DNA, RNA and Proteins
without compromising the integrity of your sample!



1.800.932.5000 | vwr.com



Prices and product details are current when published; subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks, of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s). | Visit vwr.com to view our privacy policy, trademark owners and additional disclaimers. ©2016 VWR International, LLC. All rights reserved.

0916 Lit. No. 131809W

